



Merge Eye Care PACS™

A web-based image review and management system

Vendor neutral
PACS that **brings
hundreds of
modalities** into
one viewer and
communicates to
any healthcare
IT system through
DICOM and HL7

You need a PACS that is Meaningful Use Stage 2 (MU2) certified for EHR integration so that you can quickly and easily meet imaging requirements. Your practice has a vast amount of image-intensive data from various diagnostic devices used by a range of specialists. You would like to have a system that integrates all of your patient information into a single viewer so you can accommodate higher patient volumes and improve the quality of care.

Now you can.

How We Help You

Merge Eye Care PACS, modularly certified for MU2, lets you automatically import all of the images and reports from multiple diagnostic devices into a single web-based viewer and review them anywhere, any time.

Merge Eye Care PACS

is an advanced ophthalmic image management system.

With Merge Eye Care PACS you can...



Have full image viewing and manipulation functionality in a web browser



Reduce clicks through SmartViews, an innovative way of opening multiple study types, arranging them on the screen, and using Autohistory to automatically recall



Maximize interoperability by leveraging industry standards of DICOM, HL7, and the cloud for maximum portability and compatibility of image records



Integrate directly with EHRs for viewing patient studies and reports in context



View and manipulate DICOM OCT images, eliminating the need to purchase separate review licenses from modality vendors

Merge Eye Care Solution Suite

Merge Eye Station™

Utilize your existing fundus camera to a fully digital system

Merge Eye Care PACS™

Automatically import images and diagnostic reports from diagnostic devices into a single web-based viewer

Merge Eye Care Solutions

Merge Eye Care Solutions provide a suite of products specifically for ophthalmic imaging and interoperability, including a digital imaging PACS and a modular imaging device that can capture images of both anterior and posterior segments.

Merge Eye Care Solutions provide an opportunity for practices to meet the growing need for increased patient visits, early detection and accurate tracking of disease progression.

Merge Eye Care PACS is not FDA-cleared for diagnostic use on mobile devices.

This EHR Module is 2014 Edition compliant and has been certified by an ONC-ACB in accordance with the applicable certification criteria adopted by the Secretary of the U.S. Department of Health and Human Services. This certification does not represent an endorsement by the U.S. Department of Health and Human Services or guarantee the receipt of incentive payments. Drummond Group is accredited by ANSI and approved by ONC for the ONC HIT Certification Program to certify: Complete EHR, EHR Module (all), and Certification of other types of HIT for which the Secretary has adopted certification criteria under Subpart C of 45 CFR.

SPY-5699 Rev 4.0

MERGE
An IBM Company

877.446.3743 x3 • merge.com

 @MergeHealthcare

 [linkedin.com/company/merge-healthcare](https://www.linkedin.com/company/merge-healthcare)

 [facebook.com/MergeHealthcare](https://www.facebook.com/MergeHealthcare)