

Merge Hemo™

SERVER AND WORKSTATION SECURITY PATCH TESTING

Merge Hemo 9.x, 10.x Server Security Patch Testing

Security Patch	Compatible?
Windows Server 2012 R2, SQL Server 2014 – All Microsoft Patches up until 1/1/201 have been tested and can be applied.	YES
Windows Server 2012 R2, SQL Server 2008 SP3, SQL Server 2008 SP4, Cumulative Update 12 – All Microsoft Patches up until 1/1/2018 have been tested and can be applied.	YES
Windows Server 2008 R2, SQL Server 2008 SP3, SQL Server 2008 SP4, Cumulative Update 12 – All Microsoft Patches up until 1/11/2018 have been tested and can be applied.	YES
Windows Server 2003, SQL Server 2008 SP3, SQL Server 2008 SP4, Cumulative Update 12 - All Microsoft Patches up until 7/20/2015 have been tested and can be applied, except as noted below: (Note: no longer supported by Microsoft)	YES
Intel - Network - Intel® PRO/1000 MT Desktop Adapter	NO

NOTE: For Windows Server 2003 this information also applies to versions 6.10 to 8.10 since this system is an upgrade from all prior versions.

Merge Hemo 9.x, 10.x Workstation Security Patch Testing

Security Patch	Compatible?
Windows 10 - All Microsoft Patches up until 1/11/2018 have been tested including, IE11 and can be applied.	YES
Windows 7 and Windows 8 Workstations - All Microsoft Patches up until 12/13/2017 have been tested including Win 7 SP1, IE11 and can be applied, except as noted below:	YES
nVidia - Display, Other hardware - NVIDIA Quadro 600 Published 3/9/2012	NO
Windows XP Workstations - All Microsoft Patches up until 7/15/2014 have been tested and can be applied, except as noted below (Note: no longer supported by Microsoft)	YES
Intel - Network - Intel® PRO/1000 MT Desktop Adapter	NO
nVidia - Display, Pther hardware - NVIDIA Quadro NVS 285	NO

NOTE: For Windows XP Workstations this information also applies to versions 6.10 to 8.10 since this system is an upgrade from all prior versions.

©2018 Merge Healthcare Incorporated. All rights reserved. This document and any page thereof may not be copied, distributed or otherwise reproduced or electronically shared without the advance written permission of Merge Healthcare Incorporated.