

Merge OrthoTemplates™

V. 40.1

RELEASE NOTES

Merge Healthcare
900 Walnut Ridge Drive,
Hartland, WI 53029
USA

© Copyright 2017 Merge Healthcare Incorporated, an IBM Company.

The content of this document is confidential information of Merge Healthcare Incorporated and its use and disclosure is subject to the terms of the agreement pursuant to which you obtained the software that accompanies the documentation.

Merge Healthcare® is a registered trademark of Merge Healthcare Inc.

The Merge Healthcare logo is a trademark of Merge Healthcare Inc.

All other names are trademarks or registered trademarks of their respective companies.

U.S. GOVERNMENT RESTRICTED RIGHTS:


This product is a "Commercial Item" offered with "Restricted Rights." The Government's rights to use, modify, reproduce, release, perform, display or disclose this documentation are subject to the restrictions set forth in Federal Acquisition Regulation ("FAR") 12.211 and 12.212 for civilian agencies and in DFARS 227.7202-3 for military agencies. Contractor is Merge Healthcare.

INTENDED USE:

Two and three dimensional image review, manipulation, analysis, and therapy planning capabilities that support image management display needs in the medical environment from multiple locations within and outside the hospital.

The symbols glossary is provided electronically at <http://www.merge.com/Support/Resources.aspx>.

CAUTION: U.S. federal law restricts this device to sale by, or on the order of, a physician.

 **Manufacturer's Address**
Merge Healthcare Incorporated
900 Walnut Ridge Drive
Hartland, WI 53029

For application support or to report issues with user documentation, contact Customer Support:

- In North America: call toll-free 1-877-741-5369
- Outside of North America: +31.40.299.0773
- Email: MergeSupport@us.ibm.com

Part	Date	Revision	Description
TMPL-3998	July 2017	1.0	Initial revision for version 40.1

About the Library

Implant overlays are two-dimensional templates representing prosthetic implants used in the planning of orthopaedic surgery. The digital Merge OrthoTemplates library contains orthopaedic implant overlay templates from several manufacturers. These prosthetic templates can be imported into Merge Healthcare applications to assist with surgical planning and implant selection. For a list of templates, see [“Templates Available in Version 40.1”](#). To obtain this latest library, contact [Customer Support](#).

The minimum supported version of Merge OrthoCase is version 3.6.

Merge does validation of each OrthoTemplate Library release with the latest version of OrthoCase available.

As of OrthoCase 3.7.1, multi-line text labels (as encoded by Templates manufacturers) can be displayed. Templates with multi-line text cannot be used with prior versions of OrthoCase.

Features

The following features are available with the Merge OrthoTemplates library.

Orthopaedic Templating

After you calibrate your image and select a Merge Healthcare implant overlay, the application automatically scales the implant overlay to the true image size based on the image calibration.

Orthopaedic Automatic Hip Templating

With the Merge OrthoTemplates and *ImageSnap Technology™*, Merge Healthcare introduces cutting-edge orthopaedic image processing technology. ImageSnap aims to streamline the orthopaedic workflow by automatically detecting anatomical structures. Manual tasks such as implant size selection and placement are performed automatically by the application using sophisticated matching algorithms.

Implant overlays that support the Automatic Hip Templating feature are highlighted in blue. If you select an implant family that supports this feature, the application automatically selects a hip implant size and places the hip implant overlay in a suggested location.

Templates Available in Version 40.1

The following table lists the implant overlay templates. Templates that have been deleted are indicated in black. Templates that have been updated, renamed, or added are highlighted. Their specific components and part numbers are detailed within the Merge OrthoCase application.

Products named in this document may be the trademarks or registered trademarks of their respective companies.

Manufacturer	Category	Implant Overlay Template Family
Aesculap	Total Hips	Bicontact D Plasmapore 8/10
		Bicontact E Cementless Stem
		Bicontact EH Cementless Stem
		Bicontact N Plasmapore 8/10
		Bicontact S Plasmapore 12/14
		Bipolar Cup - Inner Diameter - 22.2mm
		Bipolar Cup - Inner Diameter - 26mm
		Bipolar Cup - Inner Diameter - 28mm
		Excia 8/10 cementless
		Excia L 8/10 cementless
		Modular Head 12/14 32mm - Ceramic
		Modular Head 12/14 32mm - Metal
		Modular Head 8/10 22.2mm - Ceramic
		Modular Head 8/10 22.2mm - Metal
		Modular Head 8/10 26mm - Ceramic
		Modular Head 8/10 26mm - Metal
		Modular Head 8/10 28mm - Ceramic
		Modular Head 8/10 28mm - Metal
		Modular Head 8/10 32mm - Ceramic
		Modular Head 8/10 32mm - Metal
		Modular Head Biolox 12/14
		Modular Head COCR 12/14
		Plasmacup Insert 22.2mm Posterior Wall
		Plasmacup Insert 22.2mm Symmetrical (Polyethylene)
		Plasmacup Insert 26mm Posterior Wall
		Plasmacup Insert 26mm Symmetrical (Polyethylene)
		Plasmacup Insert 28mm Posterior Wall
		Plasmacup Insert 28mm Symmetrical (Ceramic)
		Plasmacup Insert 28mm Symmetrical (Polyethylene)
		Plasmacup Insert 32mm Posterior Wall
		Plasmacup Insert 32mm Symmetrical (Ceramic)
		Plasmacup Insert 32mm Symmetrical (Polyethylene)
		Plasmacup SC
		Trilliance Triple Tapered Polished Hip Stem
		TRJ High Offset Stem with 123° CCD
		TRJ Standard Offset Stem with 131° CCD
		Total Knees
	Columbus Femur CR/RP	
	Columbus Femur PS	
	Columbus Tibial Plateau	
	e.motion Femoral Component	
	e.motion Patella	
	e.motion Tibial Plateau	
Vega Femur PS		
Vega Patella		

Manufacturer	Category	Implant Overlay Template Family	
		Vega Tibia PS	
		Vega Tibia PS+	
		Vega Tibial Stems cemented	
		Vega Tibial Stems cementless	
	Trauma	Targon PF Long Nails	
		Targon PF Nails	
		Targon Proximal Nails	
		Targon RF Nail with Interlocking Screw (dia 6mm)	
		Targon RF Nails	
		Targon Shaft Nails	
	Arthrex	Total Knees	iBalance TKA Femoral Component
			iBalance TKA Tibial Component
		Total Shoulders	Univers(TM) Keeled Glenoids (new)
Univers(TM) Pegged Glenoids (new)			
Univers Apex Stems			
Univers CA Humeral Heads			
Univers Fracture System Head			
Univers Fracture System Stem			
Univers II Total Shoulder System Head			
Univers II Total Shoulder System Stem			
Univers REVERS Revision Stem			
Univers REVERS Stem			
Univers REVERS Glenoid Baseplate			
UniCondylar Knees		iBalance UKA Femoral Component	
		iBalance UKA Tibial Tray Component	
Biomet	Total Elbows	Comprehensive SRS-Discovery Elbow Ulna Components	
		Discovery Elbow / Humeral	
		Discovery Elbow System - Humeral Component with Bond Coating	
		Discovery Elbow System - Ulna Component with Bond Coating	
		Kudo Elbow / Humeral	
		Kudo Elbow / Ulnar	
	Total Hips	Acetabulum CMK	
		Answer Hip System	
		AR-COM All-Poly Cemented Acetabular	
		Arcos broached proximal body	
		Arcos calcar +0 proximal body	
		Arcos calcar +10 proximal body	
		Arcos calcar +20 proximal body	
		Arcos calcar proximal body	
		Arcos Cone proximal body	
		Arcos modular bowed distal femoral components with bullet tip	
		Arcos modular bowed distal femoral components with slots	
		Arcos modular bowed interlocking distal femoral components	
		Arcos modular STS distal femoral components	
		Arcos modular STS ETO distal femoral components	
		Bi-Metric Head/Neck Femoral	
		Bi-Metric Hip Systems - Collared	
		Bi-Metric Hip Systems - Collarless	
Bi-Metric Long Stem			
Bi-Metric XR Series – CDH Femoral 131.5°			
Bi-Metric XR Series – RPP Femoral 126.5°			

Manufacturer	Category	Implant Overlay Template Family
		Bi-Metric XR Series – RPP Femoral 131.5°
		Bi-Metric XR Series – SPP Femoral 131.5°
		Bio-Moore II Taper Type I
		Bio-Moore II Taper Type I - Straight Fenestrated
		Calque CMK Dysplasique Ø22.2
		Calque CMK Original Concept Ø22.2
		Calque CMK Original Concept Ø28
		CMK Acetabular Cross Shell
		Distal PMMA Positioner
		Echo Bi-Metric Collarless Femoral,FPP LAT Offset
		Echo Bi-Metric Collarless Femoral,FPP STD Offset
		Echo Bi-Metric Collarless Femoral,RPP LAT Offset
		Echo Bi-Metric Collarless Femoral,RPP STD Offset
		Exceed ABT E1 Cemented Cup System - Flanged
		Exceed ABT E1 Cemented Cup System - Non-Flanged
		G7 Acetabular Shell
		G7 Bispherical Acetabular Shells HA Coated
		G7 Bispherical Acetabular Shells Non-Plugged
		G7 Bispherical Acetabular Shells Pre-Plugged
		Generation 4 Hip System – Cemented Femoral Components
		GTS Standard Stem
		GTS Varized Stem
		Integral Porous Primary Femoral Component
		Integral Centralizer Femoral Component
		Intrigue Cemented Femoral Component
		Lateral trochanteric claws
		M²A-38 Acetabular Cup
		M²A-Magnum Metal to Metal Acetabular Cup
		M²A-Taper Radial 2 Hole
		M²A-Taper Universal Two Hole Shells
		Mallory Head Hip Program Radial Acetabular Components
		Mallory/Head Modular Bowed Distal Femoral Components
		Mallory/Head Modular Calcar Replacement 165mm
		Mallory/Head Modular Head/Neck
		Mallory/Head Modular STS Distal Femoral Component
		Mallory/Head Protrusio Cage
		Mallory/Head Radial Acetabular Shell
		Mallory/Head Revision Femoral Components
		Mallory/Head XR Femoral Component
		Mallory/Head XR Femoral Component Lateralized
		Mallory/Head X-Series Femoral Component
		Mallory-Head Hip System, Four Fin, Three Hole RingLoc Acetabular
		MaxTi Acetabular Cage
		McLaughlin +5 Ring-Loc Acetabular
		Modular Reach Porous Proximal Segment Collared Hip
		Modular Reach Porous Proximal Segment Collarless Hip
		Modular Stanmore Hip Standard Stem
		Modular Stanmore Hip Straight Stem
		Orthopedic Slavage System Total Femur Couplers
		Par 5 Acetabular Shell
		PPS +3 Liners

Manufacturer	Category	Implant Overlay Template Family
		PPS +5 Liners
		PPS Cups
		PPS Standard
		Progressive Femoral Press Fit
		Ranawat Burnstein Acetabular Cup
		Ranawat Burnstein Ringloc
		Reach Revision Hip System 200mm Bowed Femoral Implant
		Reach Revision Hip System 200mm Straight Femoral Implant
		Reach Revision Hip System 250mm Bowed Femoral Implant
		ReCap Resurfacing Head
		Recovery Protusio Cage
		Regenerex Porous Metal Cup
		REGENEREX RINGLOC+ ACETABULAR SHELL
		Ringloc Bi-Polar Acetabular Cup
		Ringloc Hemispherical Acetabular Cup - Porous
		Rx90 Cemented Hip System
		Rx90 Revision Hip System
		Sirius Cemented Stem System
		Taperloc Complete Hip System Full Length XR 123° Femoral Component Reduced Distal, Type 1 Tapers
		Taperloc Complete Hip System Full Length XR 123° Femoral Component, Type 1 Tapers
		Taperloc Complete Hip System, 12/14 Tapers, High Offset
		Taperloc Complete Hip System, 12/14 Tapers, Standard Offset
		Taperloc Complete Hip System, Bone Master, 12/14 Tapers, High Offset
		Taperloc Complete Hip System, Bone Master, 12/14 Tapers, Standard Offset
		Taperloc complete Hip System, Bone Master, Reduced Distal 12/14 Tapers, High Offset Femoral Component
		Taperloc complete Hip System, Bone Master, Reduced Distal 12/14 Tapers, Standard Offset Femoral Component
		Taperloc Complete Hip System, Bone Master, Reduced Distal, Type 1 Tapers, High Offset
		Taperloc Complete Hip System, Bone Master, Reduced Distal, Type 1 Tapers, Standard Offset
		Taperloc Complete Hip System, Bone Master, Type 1 Tapers, High Offset
		Taperloc Complete Hip System, Bone Master, Type 1 Tapers, Standard Offset
		Taperloc Complete Hip System, Type 1 Tapers, High Offset
		Taperloc Complete Hip System, Type 1 Tapers, Reduced Distal, High Offset
		Taperloc Complete Hip System, Type 1 Tapers, Reduced Distal, Standard Offset
		Taperloc Complete Hip System, Type 1 Tapers, Standard Offset
		Taperloc Complete Microplasty Hip System High Offset Femoral Component
		Taperloc Complete Microplasty Hip System Microplasty XR 123° Femoral Component Reduced Distal, Type 1 Tapers
		Taperloc Complete Microplasty Hip System Microplasty XR 123° Femoral Component, Type 1 Tapers
		Taperloc Complete Microplasty Hip System Standard Offset
		Universal Ringloc Acetabular System
		Vision Acetabular Components
	Total Knees	AGC 3000 PS Knee System - Interlok
		AGC 3000 PS Knee System - Porous
		AGC Total Knee System Universal Femoral
		AGC Total Knee System Universal Femoral - Interlok
		AGC Total Knee System Universal Femoral - Porous
		AGC Total Knee System, Anatomic Femoral - Interlok
		AGC Total Knee System, Anatomic Femoral - Porous
		AGC Tradition Series High Post Posterior Stabilized Component Interlok Molded Tibial
		AGC Tradition Series High Post Posterior Stabilized Components, Interlok Universal Femoral
		AGC Universal PS Knee System
		Ascent Knee System - Constrained Bearing

Manufacturer	Category	Implant Overlay Template Family
		Ascent Knee System - PS Bearing
		Ascent Knee System Femoral Revision Femoral Component
		Ascent Knee System Femoral Revision Stem Extension
		Ascent Knee System Primary Curved Bearing
		Ascent Knee System Primary Lipped Bearing
		Ascent Knee System Tibial Revision
		Biomet Primary Tibial Tray
		Biomet Vanguard Knee Posterior Stabilized Bearing
		Biomet Vanguard Knee Primary Bearing
		Biomet Vanguard Knee Primary Lipped Bearing
		Biomet Vanguard Knee PS Plus Bearing
		Compress Adapter
		Compress Centering Sleeve
		Compress Transverse Pins
		Finned Stem
		Finned Tibial Tray
		I-Beam Stem
		I-Beam Tibial Tray
		Maxim Knee Constrained Bearing
		Maxim Knee Primary Femoral – Interlock Component
		Maxim Knee Primary Femoral – Porous Component
		Maxim Knee Primary Tibial Bearing
		Maxim Knee Revision Stem
		Maxim P.S Femoral Component
		Maxim Revision Femoral Component
		Microplasty Tibial Tray
		Offset Tibial Adaptor
		Offset Tibial Tray
		Offset Tray System
		OSS 10/20mm Tibial Block Augments
		OSS 3cm Non-Modular Proximal Tibia
		OSS 5cm Non-Modular Proximal Tibia
		OSS 7cm Non-Modular Proximal Tibia
		OSS 9cm Modular Proximal Tibia
		OSS Anterior Flange Augment
		OSS Bowed Collared Stem
		OSS Cemented Bowed IM Stem
		OSS Cemented Bowed Stem 150mm
		OSS Cemented Bowed Stem 225mm
		OSS Cemented Straight Stem 150mm
		OSS Cemented Straight Stem 225mm
		OSS Cemented Straight Stem 90mm
		OSS Compress Knee
		OSS Diaphyseal Segment
		OSS Elliptical Diaphyseal Segment
		OSS Elliptical Segmental Femoral
		OSS Femoral Sleeve Augment
		OSS Knee Femoral Component
		OSS Modular Tibial Base Plate
		OSS Modular Tibial Base Plate Sleeve Augments
		OSS Non-Modular Tibial Base Plate - Long

Manufacturer	Category	Implant Overlay Template Family
		OSS Non-Modular Tibial Base Plate - Short
		OSS Porous Bowed Stem 150mm
		OSS Porous Bowed Stem 225mm
		OSS Porous Coated Bowed IM Stem
		OSS Porous Straight Stem 150mm
		OSS Porous Straight Stem 90mm
		OSS Proximal Femoral
		OSS Proximal Tibial Sleeve
		OSS Resurfacing Femoral
		OSS RS Elliptical Segmental Femoral
		OSS RS Proximal Tibial Sleeve
		OSS RS Resurfacing Femoral
		OSS RS Segmental Femoral
		OSS Segmental Adapter
		OSS Segmental Femoral
		OSS Straight Collared Stem
		Oxford Partial Knee Twin-Peg
		Oxford Unicompartmental Knee System
		PS Tibial Bearing Revision Component
		PSC Tibial Bearing Revision Component
		Repicci II Knee Tibial Component
		Short Compress 30 mm Anti-Rotation Spindle
		Short Compress 38 mm Anti-Rotation Spindle
		Short Compress 44 mm Anti-Rotation Spindle
		Short Compress Elliptical Anti-Rotation Spindle
		Standard Compress 30 mm Anti-Rotation Spindle
		Standard Compress 38 mm Anti-Rotation Spindle
		Standard Compress 44 mm Anti-Rotation Spindle
		Standard Compress Elliptical Anti-Rotation Spindle
		Stemed Tibial Tray
		Vanguard Complete Knee System Primary PS Femoral Component
		Vanguard Knee Open Box PS Femoral
		Vanguard Knee Primary CR Femoral - Interlock
		Vanguard Knee Primary CR Femoral - Porous
		Vanguard Knee SSK INTLK Femoral
		Vanguard Knee SSK PS Tibial Bearing – 59mm Revision Tibial Component
		Vanguard Knee SSK PS Tibial Bearing – 63mm/67mm Revision Tibial Component
		Vanguard Knee SSK PS Tibial Bearing – 63mm/67mm Revision Tibial Component – Small Post
		Vanguard Knee SSK PS Tibial Bearing – 71mm/75mm Revision Tibial Component
		Vanguard Knee SSK PS Tibial Bearing – 71mm/75mm Revision Tibial Component – Small Post
		Vanguard Knee SSK PS Tibial Bearing – 79mm/83mm Revision Tibial Component
		Vanguard Knee SSK PSC Tibial Bearing – 59mm Revision Tibial Component
		Vanguard Knee SSK PSC Tibial Bearing – 63mm/67mm Revision Tibial Component
		Vanguard Knee SSK PSC Tibial Bearing – 63mm/67mm Revision Tibial Component – Small Post
		Vanguard Knee SSK PSC Tibial Bearing – 71mm/75mm Revision Tibial Component
		Vanguard Knee SSK PSC Tibial Bearing – 71mm/75mm Revision Tibial Component – Small Post
		Vanguard Knee SSK PSC Tibial Bearing – 79mm/83mm Revision Tibial Component
		Vanguard RP Hi-Flex Femoral
		Vanguard RP Tibial Tray with Vanguard RP Hi-Flex Bearing
		VANGUARD SSK 360 Femoral Component
		VANGUARD SSK 360 Revision Knee System - Bowed Splined

Manufacturer	Category	Implant Overlay Template Family
		VANGUARD SSK 360 Revision Knee System - Grit Blast Bowed
		VANGUARD SSK 360 Revision Knee System - Grit Blast Straight
		VANGUARD SSK 360 Revision Knee System - Offset Adapters
		VANGUARD SSK 360 Revision Knee System - Smooth Bowed
		VANGUARD SSK 360 Revision Knee System - Smooth Straight
		VANGUARD SSK 360 Revision Knee System - Straight Splined
		VANGUARD SSK 360 Revision Knee System - Tibial Component
		VANGUARD SSK 360 Revision Knee System - Tibial Keel Augment
		Total Shoulders
	ArComXL Humeral Bearing With Humeral Tray	
	Bi-Angular Humeral Stem	
	Bi-Angular Humeral Stem - Absolute Bi-Polar Head	
	Bio-Modular Humeral Stem - Absolute Bi-Polar Head	
	Bio-Modular Shoulder - Humeral Head Components	
	Bio-Modular Shoulder - Porous Coated Glenoids - Porous Coated Tray & Arcom liner	
	Bio-Modular Shoulder - Primary Humeral Stem	
	Bio-Modular Shoulder - Revision Humeral Stem	
	Comprehensive Reverse Shoulder System - Glenosphere Baseplate W/ Taper Adaptor	
	Comprehensive Reverse Shoulder System - Versa-Dial Glenosphere	
	Comprehensive Revision Stems	
	Comprehensive Shoulder System – Humeral Fracture Stem	
	Comprehensive Shoulder System – Humeral Primary Stem	
	Comprehensive SRS - Bio-Modular Shoulder Humeral Heads	
	Comprehensive SRS - Comprehensive Reverse Shoulder System Heads	
	Comprehensive SRS - Distal Bodies Right and Left	
	Comprehensive SRS - Extended Articular Surface Humeral Heads	
	Comprehensive SRS - Intercalary Segments	
	Comprehensive SRS - Large Revision Proximal Body	
	Comprehensive SRS - Small Revision Proximal Body	
	Comprehensive SRS - Stems	
	Comprehensive SRS - Super Extended Articulation Surface Head	
	Comprehensive SRS - Tumor Proximal Body	
	Comprehensive SRS - Versa-Dial Head	
	Mini Compress/Mosaic Taper Adapter - Mini	
	Mini Compress/Mosaic Taper Adapter - Standard	
	Mini Compress/SRS Taper Adapter - Mini - 30°	
	Mini-Compress Mosaic adapter Long Anchor Plug	
	Mini-Compress Mosaic adapter Short Anchor Plug	
	Mosaic Humeral Replacement - IM Stem System	
	Mosaic Humeral Replacement System - Humeral Head	
	Mosaic Humeral Replacement System - Proximal Body And Segment	
	Mosaic Humeral Replacement System - Soft Tissue Attachment Sleeve	
Short Mini Compress 25 mm Spindle		
Short Mini Compress 30 mm Anti-Rotation Spindle		
Standard Compress/SRS Taper Adapter Standard - 30°		
Standard Mini Compress 25 mm Spindle		
Standard Mini Compress 30 mm Anti-Rotation Spindle		
Versa-Dial Humeral Head		
Trauma	DVR Anatomic Plate	
	AFFIXUS Hip Fracture 1.5M Nail System - Anti-Rotation Screws	
	AFFIXUS Hip Fracture 1.5M Nail System - Distal Screws	

Manufacturer	Category	Implant Overlay Template Family
		AFFIXUS Hip Fracture 1.5M Nail System - End Cap Flush
		AFFIXUS Hip Fracture 1.5M Nail System - In Cap Flush
		AFFIXUS Hip Fracture 1.5M Nail System - In Cap Flush Impinging
		AFFIXUS Hip Fracture 1.5M Nail System - Lag Screws
		AFFIXUS Hip Fracture 1.5M Nail System - Long Nail 125°
		AFFIXUS Hip Fracture 1.5M Nail System - Short Nail 125°
		AFFIXUS Hip Fracture 1.5M Nail System - Short Nail 130°
		UniCondylar Knees
	Oxford Partial Knee High Flex Femoral	
	Oxford Partial Knee Standard Femoral	
	Oxford Partial Knee Tibial Tray	
	Vanguard M Series Tibial Bearings	
	BIOTECHNI	Total Hips
HIPFORTRESS-ND cementless stem 130° collarless		
Snapfit Bipolar Cup		
Snapfit Modular Head		
Corin	Total Hips	CTIII Cementless Stem
		CoCr Modular Head Cone 11/13
		Cormet Acetabular Cup
		Cormet Head
		Metafix Cementless Lateralised Stem 135°
		Metafix Cementless Standard Stem 125°
		Metafix Cementless Standard Stem 135°
		MiniHip Cementless Femoral Stem 11/13 Taper
		TRIFIT LATERALISED TAPERED STEM
		TRIFIT STANDARD TAPERED STEM
DePuy Johnson & Johnson	Total Hips	AML A Plus Hip Stem
		AML Hip
		Arthopor II Oblong Porous Coated Acetabular Cup 15mm
		Arthopor II Oblong Porous Coated Acetabular Cup 25mm
		ASR Acetabular Cup
		ASR UNI Femoral
		Charnley Flanged Cup
		Charnley LPW Cup
		Charnley Modular Stem
		Charnley Ogee Cup
		Charnley Standard Cup
		Charnley Stem
		CORAIL AMT COXA VARA
		CORAIL AMT DYSPLASIA COLLARED
		CORAIL AMT DYSPLASIA COLLARLESS
		CORAIL AMT HIGH OFFSET
		CORAIL AMT SHORT NECK 135°, COLLARED
		CORAIL AMT SHORT NECK 135°, COLLARLESS
		CORAIL AMT STANDARD 125°, COLLARED
		CORAIL AMT STANDARD 125°, COLLARLESS
		CORAIL AMT STANDARD COLLARED
		CORAIL AMT STANDARD COLLARLESS
		Corail Stem – Collared, Coated
		C-Stem AMT
		C-Stem Asian Stem

Manufacturer	Category	Implant Overlay Template Family
		C-Stem Total Hip System
		Duraloc Acetabular Cup System
		Duraloc Option Cup
		Elite Plus CoCr Head
		Elite Plus Flanged Cup
		Elite Plus Long Posterior Wall Cup
		Elite Plus Ogee Acetabular Cup
		Elite Plus Ogee Cup
		Elite Plus Retentive Bore Flanged Cup
		Elite Plus Stainless Steel Head
		Elite Plus Stem
		Endurance Total Hip System – 28mm
		Euroform Collared Cemented Stem
		Euroform Collared Cementless Stem
		Euroform Collarless Cemented Stem
		Euroform Collarless Cementless Stem
		European Future Hip
		G2 Hip - Uncemented
		J-Fx Bipolar Cup-22mm
		J-Fx Bipolar Cup-28mm
		J-Loc Cup – 10 Degree
		J-Loc Cup – 1200 Series
		J-Loc Cup – 300 Series
		J-Loc Cup – Neutral
		Kar Revision Corail Stem
		Modular Cathcart Unipolar
		Pinnacle Acetabular Cup
		Pinnacle Acetabular Cup System
		Prodigy Hip System
		Protrusio Cage
		Reclaim Revision Hip System, Conical Proximal Body
		Reclaim Revision Hip System, Tapered Distal Stems
		Reef Collared/Collarless Trochanteric Component
		Reef Lateral Wing
		Reef Stem
		Replica A Stem
		Self-Centering Bipolar
		Solution Acetabular Cup System
		Solution Bowed Stem
		Solution Calcar Stem
		Solution Straight Stem
		S-ROM Asian Hip Stem
		S-ROM Asian Hip System
		S-ROM Hip System – Femoral Stems Lateral
		S-ROM Hip System – Femoral Stems Offset Lateral
		S-ROM Hip System – Femoral Stems Standard
		S-ROM Hip System – Proximal Sleeves
		S-ROM Total Hip System
		Summit Basic Cemented Hip Stem - Without Measurements (renamed, was Summit Basic Cemented Hip Stem)
		Summit Basic Cemented Hip Stem - With Measurements (new)
		Summit Basic Press-Fit Hip Stem

Manufacturer	Category	Implant Overlay Template Family
		Summit Cemented Femoral Stem
		Summit Tapered Hip System
		Total ASR Femoral
		TRI LOCK Bone Preservation Stem
		Tri-Lock Total Hip System
		Ultima Cemented Stem
		Ultima Press Fit Stem
		Ultima Unipolar Head
		Vision 2000 Hip System
		Total Knees
	ATTUNE CEMENTED TIBIAL BASE ROTATING PLATFORM	
	ATTUNE CR CEMENTED FEMORAL COMPONENT	
	ATTUNE PS CEMENTED FEMORAL COMPONENT	
	LCS Complete Mobile Bearing Knee System	
	LCS Complete RP Insert	
	LCS PCR Tibial Porous Tray	
	LCS PCR Tibial Textured Tray	
	M.B.T REVISION CEMENTED TIBIAL TRAY	
	M.B.T REVISION TIBIAL SLEEVE CEMENTED	
	M.B.T REVISION TIBIAL SLEEVE POROCOAT	
	MBT Duofix Tibial Tray	
	MBT Tray System – Keel/Non-keel	
	PCR Bearing	
	PFC Sigma Adapter	
	PFC Sigma Femoral Adapter 5°	
	PFC Sigma Femoral Adapter 7°	
	PFC Sigma Knee System	
	PFC Sigma Knee System - Keel	
	PFC Sigma Knee System – Modular Plus	
	PFC Sigma RPF	
	PFC Sigma TC3 Femoral Augmentation	
	PFC Sigma Tibial Cemented Stem	
	PFC Sigma Tibial Stem Press Fit 115mm	
	PFC Sigma Tibial Stem Press Fit 150mm	
	PFC Sigma Tibial Stem Press Fit 75mm	
	PFC Tappered Cemented Stem	
	RP Bearing	
	Sigma CR 150 High Flex Knee System	
	Sigma PS 150 Femoral Component	
	S-ROM Noiles Femoral Component & Tibial Plateau	
	S-ROM Noiles Femoral Stem Extensions	
	S-ROM Noiles Modular Femoral Sleeve	
S-ROM Noiles Modular Tibial Base Component		
S-ROM Noiles Modular Tibial Sleeves		
S-ROM Noiles Tibial Cap		
S-ROM Noiles Tibial Stem Extensions		
S-ROM Noiles Tibial Tray		
Total Shoulders	CTA Head	
	Delta CTA Glenosphere	
	Delta CTA Humeral Epiphysis	
	Delta CTA Humeral Inox Heads	

Manufacturer	Category	Implant Overlay Template Family	
		Delta CTA Humeral Spacer	
		Delta CTA Reverse Shoulder System – Cemented Humeral Diaphysis	
		Delta CTA Reverse Shoulder System – Cementless Humeral Diaphysis	
		Delta CTA Reverse Shoulder System – Humeral Cup	
		Delta CTA Standard Metaglene	
		Eccentric Glenosphere	
		Global Advantage Eccentric Head	
		Global Advantage Long Humeral Stem	
		Global Advantage Standard Head	
		Global Advantage Standard Humeral Stem	
		Global AP Eccentric Heads – Humeral Head	
		Global AP Glenoids	
		Global AP Humeral Stem	
		Global AP Standard Heads – Humeral Head	
		Global Cap – Humeral Head	
		Global FX Shoulder Long	
		Global FX Shoulder Standard	
		Global Icon - Humeral Anchor	
		Global Icon - Humeral Head	
		Global Unite Anatomic Body	
		Global Unite FX Eccentric Suture Collar	
		Global Unite FX Proximal Body	
		Global Unite FX Standard Suture Collar	
		Global Unite Humeral Head Eccentric	
		Global Unite Humeral Head Standard	
		Global Unite Humeral Stem Long	
		Global Unite Humeral Stem Std	
		High Mobility Humeral PE Cup	
		Modular Humeral Stem	
		Monobloc Cemented	
		Retentive Humeral PE Cup	
		Standard Glenosphere	
		Standard Humeral PE Cup	
		Trauma	3.5mm Solid Cortical Screw, Non Sterile
			3.5mm Solid Cortical Screw, Sterile
			4.8mm Cancellous Screw, Non Sterile
			4.8mm Cancellous Screw, Sterile
			ACE DVR Anatomic Plate
			Proximal Nail
			Universal Nail
			VersaNail Femoral Troch Entry Nail
			VersaNail Femoral Universal Nail
VersaNail Tibial Nail			
UniCondylar Knees	Preservation Knee – Mobile Bearing		
	Preservation Knee – Tibial Tray		
DJO Surgical	Total Hips	FMP ACETABULAR SYSTEM - HEMI	
		FMP ACETABULAR SYSTEM - RIM FLARE	
		FOUNDATION TOTAL HIP SYSTEM	
		Linear Hip - LATERAL OFFSET	
		Linear Hip - STANDARD OFFSET	
		TaperFill Stem, Standard Offset	

Manufacturer	Category	Implant Overlay Template Family	
		TaperFill Stem, Lateral Offset	
		Total Hip - Revelation	
		Total Hip - Revelation v2	
		Total Hip - Unipolar/Bipolar Hip	
	Total Knees	3DKnee Total Knee System - Femoral Components	
		3DKnee Total Knee System - Tibial Components	
		FOUNDATION KNEE PS - FEMORAL COMPONENTS	
		FOUNDATION KNEE PS - STEMMED AND RESURFACING TIBIAL COMPONENTS	
		PRIMARY FOUNDATION KNEE - FEMORAL COMPONENTS	
		PRIMARY FOUNDATION KNEE - PRIMARY ALL POLY TIBIAL COMPONENTS	
		PRIMARY FOUNDATION KNEE - STEMMED AND RESURFACING TIBIAL COMPONENTS	
	Total Shoulders	Turon All-Poly Keeled Glenoid	
		Turon All-Poly Pegged Glenoid	
		Turon Offset Humeral Head	
		Turon Primary Humeral Stem	
		Turon Revision Humeral Stem	
		Turon Standard Humeral Head	
	Exactech	Total Hips	AcuMatch A-Series Acetabular Shell
			AcuMatch C-Series - Cemented Femoral Stems (11/13 taper), Extended Offset
AcuMatch C-Series - Cemented Femoral Stems (11/13 taper), Standard Offset			
AcuMatch M-Series Metaphyseal Calcar			
AcuMatch M-Series Metaphyseal Standard			
AcuMatch M-Series Neck			
AcuMatch M-Series Stem Curved			
AcuMatch M-Series Stem Straight			
BIPOLAR HEAD			
Novation CFS Cemented Stem			
Novation CFS Press-Fit			
Novation Crown Cup			
Novation Element, Collared, Standard Offset			
Novation Element, Collarless, Extended Offset			
Novation Element, Collarless, Reduced Offset (Extended)			
Novation Element, Collarless, Reduced Offset (Standard)			
Novation Element, Collarless, Standard Offset			
Novation Press Fit Splined Stem, Extended Offset			
Novation Press Fit Splined Stem, Standard Offset			
Novation Press Fit Tapered Stem Extended Offset			
Novation Press Fit Tapered Stem Standard Offset			
Z-Stem Lateral			
Z-Stem Revision			
Z-Stem Standard			
Total Knees		Asymmetric CR and PS	
		Constrained Condylar	
		Optetrak Logic Knee System Femoral Components	
		Optetrak Logic Knee System Tibial Tray	
		Stem Extensions	
		Symmetric	
		Tibial Finned Trays	
		Tibial Trapezoidal Trays	
		Unicondylar Knee	
Total Shoulders	38mm Reverse Humeral Components Liners +0mm		

Manufacturer	Category	Implant Overlay Template Family	
		38mm Reverse Humeral Components Liners +2.5mm	
		42mm Reverse Humeral Components Liners +0mm	
		42mm Reverse Humeral Components Liners +2.5mm	
		46mm Reverse Humeral Components Liners +0mm	
		46mm Reverse Humeral Components Liners +2.5mm	
		Cemented Humeral Fracture Stem	
		Cemented Humeral Long Stem	
		Cemented Humeral Stem	
		Cemented Platform Fracture Stem	
		Glenoids Keeled	
		Glenoids Pegged	
		Humeral Heads	
		Pressfit Humeral Stem	
		Reverse Fracture Humeral Components Liners +0	
		Reverse Fracture Humeral Components Liners +2.5	
		Reverse Glenoid Plate/Glenospheres	
		Reverse Humeral Components 15mm Tray Liners +0	
		Reverse Humeral Components 15mm Tray Liners +2.5	
Homs	Trauma	Acromioclavicular Plate	
		Clavicle Plate	
		Distal Fibula Plate	
		IPT Distal Screw	
		IPT Extra Screw	
		IPT Lag Screw	
		IPT Nail System Long Nail 125°	
		IPT Nail System Long Nail 130°	
		IPT Nail System Nail 125°	
		IPT Nail System Nail 130°	
		Reconstruction Plate	
		SBFix	
		JMMDM	Total Hips
12/14 Ovation Stem, STD			
B2 Bipolar Cup			
Ecompass Cemented Hip System			
Escalade Acetabular Cup			
Hemisphere Modular Cup			
Legend All Poly Cup			
Ovation Hip Stem STD and EXT			
Ovation Tribute Stem, EXT			
Ovation Tribute Stem, STD			
Pivot Bipolar Cup			
Primaloc Cemented Stem			
Primaloc Cementless Stem, Collared Stem			
Primaloc Cementless Stem, Collarless Stem			
Total Knees	BKS modular Cemented Stems		
	BKS Modular Femoral Component		
	BKS modular Fluted Stems		
	BKS Modular Tibial Neutral		
	BKS Modular Tibial Sloped		
	BKS Primary Femoral Component		
	BKS Primary Tibial Tray		

Manufacturer	Category	Implant Overlay Template Family
	Trauma	BKS TriMax PS Femoral Component
		Compression Hip Screw System - 4.5mm Cortical Screw
		Compression Hip Screw System - 7.0mm Cannulated Cancellous Screw
		Compression Hip Screw System - Lag Screw
		Compression Hip Screw System - Side Plate
		Compression Hip Screw System - Torochanteric plate
		Distal Fibra Plate
		Distal Radius Plate
		Distal Tibia Plate
		Locking Straight Plate
		Proximal Humeral Plate
		Reconstruction Straight Plate
		KYOCERA
TNK TIBIAL ANKLE		
Total Hips	7 Stem	
	910DTP Stem	
	910DTP Stem H/O	
	910FS Socket CP	
	910MX FS Socket CP	
	910PERFIX C Collared	
	910PERFIX HA Collared	
	910PERFIX HA Collarless	
	910PERFIX HA Collarless HO Stem	
	910PERFIX HA Collarless P+HO Stem	
	910YU Collared	
	910YU Collared H/O	
	MOD-Centaur Stem	
	MOD-Centaur Proximal Femur	
	AHFIX Q Stem NAR	
	AHFIX Q Stem NAR Collared	
	AHFIX Q Stem STD	
	AHFIX Q3 Shell with MX Liner CP - 22	
	AHFIX Q3 Shell with MX Liner CP - 26	
	AHFIX Q3 Shell with MX Liner CP - 28	
	AHFIX Q3 Shell with STD Liner CP - 22	
	AHFIX Q3 Shell with STD Liner CP - 26	
	AHFIX Q3 Shell with STD Liner CP - 28	
	AMS HA Shell 910 MX Liner	
	AMS HA Shell 910 STD Liner	
	AMS HA Shell with 910MX Liner CP - 22	
	AMS HA Shell with 910MX Liner CP - 26	
	AMS HA Shell with 910MX Liner CP - 28	
	AMS HA Shell with 910MX Liner CP - 32	
	AMS HA Shell with 910STD Liner CP - 26	
	AMS HA Shell with 910STD Liner CP - 28	
	AMS HA Shell with 910STD Liner CP - 32	
	Aquala AMS 910 MX Liner	
	Aquala AMS 910 STD Liner	
	Aquala Q3 MX Liner	
	Aquala Q3 STD Liner	
	Aquala Q5LP MX Liner	

Manufacturer	Category	Implant Overlay Template Family
		Aquala Q5LP STD Liner
		CHANGEA Stem
		DUALCUP (new)
		Elance HO
		Elance Stem STD
		INHERITOR HO
		INHERITOR STD
		INITIA HO (new)
		INITIA STD (new)
		J-Taper Stem HO
		J-Taper Stem STD
		K-MAX ABC HIP Stem Collared
		K-MAX ABC HIP Stem Collarless
		K-MAX CLHO Cup Z
		K-MAX CLHP Cup
		K-MAX DUAL CLDCS CUP
		K-MAX KT Plate S
		K-MAX S-LOCK Stem
		K-MAX SS HIP Collared
		K-MAX SS HIP Collarless
		K-MAX SS HIP Taper
		PerFix HA Collared
		PerFix HA Collarless
		PHS Outer Head 22mm
		PHS Outer Head 26mm
		SC Stem POLISHED CDH
		SC Stem POLISHED STD
		SQRUM HO
		SQRUM STD
		Superio ST Collared
		Superio ST Collarless
		Total Hip - HD Stems
		Type1 Long Stem Collarless
		TYPE6 LONG Stem
		TYPE6 Stem
		YU Collared
		YU for 26mm Head
		YU Long
	Total Knees	ACTIYAS T-TRAY ACTIYAS T-TRAY Extension 40 ACTIYAS T-TRAY Extension 75 ACTIYAS ZF Bi-Surface KU4+ Bi-Surface5 PS Bi-Surface5 T-CS Bi-Surface5 T-CS Extension 40 Bi-Surface5 T-CS Extension 75 Bi-Surface5 T-PS Bi-Surface5 T-PS Extension 40 Bi-Surface5 T-PS Extension 75 EMK CR Femoral

Manufacturer	Category	Implant Overlay Template Family	
		EMK CR Tibial	
		EMK PS Femoral	
		EMK PS Tibial	
		KU3 T-Keel NH M/TI	
		KU3 T-M/TI	
		KU4+ ZF	
		LFA CR Femoral	
		LFA CR Tibial	
		LFA T-M/TI / CR	
		LFA T-M/TI / PS	
		LFA ZF/MACRO / CR	
		LFA ZF/MACRO / PS	
		Physio-Knee F-M	
		Physio-Knee T-M/TI	
		Unicompartmental Knee Arthroplasty - Femoral Component	
		Unicompartmental Knee Arthroplasty - Tibial Component	
		Total Shoulders	PSS Ball
		PSS Glenoid	
		PSS: Plus Medium Stem	
		PSS: Short Stem	
	Trauma	AA Hip screw	
		OSR	
	Tumor	KLS DF	
		KLS Extension	
		KLS F Connector	
		KLS PF	
		KLS Stem-ARC	
		KLS Stem-ARC CUR	
		KLS Stem-CUR	
		KLS T	
	Lima Corporate	Total Hips	C2-LAT Stem
			C2-REV Stem
			C2-STD Stem
			CL-Revision Stem
		DELTA-PF Cup	
		DELTA-REVISION Cup	
		DELTA-REVISION-TT Cup	
		DELTA-ST_C Cup	
		DELTA-TT Cup	
		Friendly Hip Stem	
		H-MAX S Lateralizing Stem (new)	
		H-MAX S Standard Stem (new)	
		Lock Bipolar AISI 316/L Heads	
		Lock Bipolar CoCrMo Heads	
		Master SL Lateralizing Stem - Poro Ti Coating	
		Master SL Standard Stem - Poro Ti Coating	
		MINIMA S LATERALIZING Stem, TAPER 12/14	
		MINIMA S STANDARD Stem, TAPER 12/14	
		Modular Neck	
		Modular Stem	
		Modulus-R femoral stem	

Manufacturer	Category	Implant Overlay Template Family
		PROTRUTED Acetabular Cup
		REVISION DELTA-ONE-TT Cup
		REVISION Hip Cementless Modular stem
		Revision Lateralized Neck (new)
		Revision Modular Stem (new)
		Revision Neck (new)
		STANDARD Acetabular Cup
	Total Shoulders	SMR Cemented Glenoid
		SMR Cemented Stem
		SMR Cementless Finned Stem
		SMR Cementless Glenoid L1 Plus UHMWPE Liner
		SMR CTA Humeral Head
		SMR Eccentric Glenosphere Plus Glenoid
		SMR Finned Humeral Body
		SMR Glenosphere Plus Glenoid
		SMR Humeral Head
		SMR Reverse Humeral Body Plus Retentive Cross-linked UHMWPE Liner
		SMR Reverse Humeral Body Plus Retentive UHMWPE Liner
		SMR Reverse Humeral Body Plus Reverse Cross-linked UHMWPE Liner
		SMR Reverse Humeral Body Plus Reverse UHMWPE Liner
SMR Trauma Humeral Body		
Trauma	Acetabular Plate for Cemented Cups-Gabbie per Cotili da Revisione	
LINK	Total Hips	MP Neck Segment without collar (new)
		MP Neck Segment XXL without collar (new)
		MP Reconstruction Collar
		MP Reconstruction Neck Segment
		MP Reconstruction Stem
Mathys	BiCondylar Knees	balanSys BiCondylar Femur
		balanSys Tibia Plateau
	Total Hips	Bipolar Head CoCrMo 22.2
		Bipolar Head CoCrMo 28
		CBC Stem - Standard
		CBH Stem - Lateral
		CBH Stem - Standard
		CBK Revision Stem
		CCA Stem - Lateral
		CCA Stem - Standard
		CCE Roof Reinforcement Ring
		Centris Displasy Stem
		Centris Long Stem
		Centris Standard Stem
		expanSys MoM
		expanSys PE Inlay Shoulder
		expanSys PE Inlay Standard
		Full Profile CCB Cup
		Full Profile RM Cup
		GSS Compulsory Cemented Stem
		GSS-CL Cementless Stem
		GSS-CO Cementless Stem
		Hemiprosthesis Bipolar Head
		Hemiprosthesis Ceramic Head

Manufacturer	Category	Implant Overlay Template Family
		Hemiprosthesis Stainless Steel Head
		Low Profile CCB Cup
		Modular Revision Head/Neck
		Modular Revision Stem
		Optimys Stem Uncemented Lateral
		Optimys Stem Uncemented Standard
		RM Bevelled Cup
		RM MoM
		RM Pressfit
		RM Pressfit MoM
		RM Pressfit vitamys
		RM Revision
		seleXys CoC PC
		seleXys CoC TH
		seleXys CoC TH+
		seleXys CoC TPS
		seleXys CoC TPS-R
		seleXys MoM PC
		seleXys MoM TH
		seleXys MoM TH+
		seleXys MoM TPS
		seleXys MoM TPS-R
		seleXys Poly 0° PC
		seleXys Poly 0° TH
		seleXys Poly 0° TH+
		seleXys Poly 0° TPS
		seleXys Poly 0° TPS-R
		seleXys Poly 10° PC
		seleXys Poly 10° TH
		seleXys Poly 10° TH+
		seleXys Poly 10° TPS
		seleXys Poly 10° TPS-R
		Standard Stem
		twinSys Cemented Lateral
		twinSys Cemented Standard
		twinSys Uncemented Lateral
		twinSys Uncemented Long Stem
		twinSys Uncemented Standard
		twinSys XS
	Total Knees	balanSys Femur CR
		balanSys Femur PS
		balanSys Tibia CR
		balanSys Tibia PS CR
		balanSys Tibia PS
	Total Shoulders	Affinis Cemented Stem
		Affinis CoCr Head
		Affinis Fracture Connector
		Affinis Fracture Inverse
		Affinis Fracture Stem
		Affinis Head Adaptor
		Affinis Inverse Cemented Stem

Manufacturer	Category	Implant Overlay Template Family
		Affinis Inverse Inlays
		Affinis Inverse Metaglene-Glenosphere
		Affinis Inverse Revision Metaglene-Glenosphere
		Affinis Inverse Uncemented Stem
		Affinis Short
		Affinis Short Stem
		Affinis Uncemented Stem
		Articula Connector
		Articula Inverse
		Articula Stem
		UniCondylar Knees
Matortho	Total Hips	ADEPT Extra Fixation Metal On Metal Cup
		ADEPT Metal On Metal Cup
		ADEPT Metal on Metal Resurfacing Head
	Total Knees	Saiph Knee Femoral Component
		Saiph Knee Tibial Component
Medacta	Total Hips	AMISem C
		AMISem H
		AMIS-K Cemented Stem (new)
		Femoral Head 12/14 CoCr (new)
		Femoral Head 12/14 MectaCer BioloX Delta (new)
		Femoral Head 12/14 MectaCer BioloX Option System (new)
		MasterLoc Mectagrip Coated (updated)
		Mpact Multiholes (new)
		Mpact Two Holes (new)
		Quadra C
		Quadra H
		Quadra S
		Quadra SN C
	Quadra SN H	
	Total Knees	GMK-Primary PS Femoral Component
		GMK-Primary STD Femoral Component
		GMK-Primary Tibial Tray Fixed
		GMK-SPHERE Primary Femoral Component
		GMK-SPHERE Tibial Tray Fixed
Microport	Total Hips	PROFEMUR(R) AM
Mizuho	Total Hips	ARC HA System CH 28 MX
		ARC HA System CH 28 STD
		ARC HA System CH 32 MX
		ARC HA System CH 32 STD
		ARC HA System MH 28 MX
		ARC HA System MH 28 STD
		ARC HA System MH 32 MX
		ARC HA System MH 32 STD
		AW-GC Cup 22EV
		AW-GC Cup 26EV
		AW-GC Cup 26FL
		Cement Socket 22
		Cement Socket 26
		DK Bipolar
		MaximLock

Manufacturer	Category	Implant Overlay Template Family
		MX-10 Hip System
		MX-200 Hip System
		MX-2000 Hip System
		MX-2000 Locking Screw
	Trauma (new)	HYBRIX Plate (new)
Nemoto	Total Hips	Y Shell
		Y Stem
NexMed International	Total Hips (new)	Bipolar Cup (new)
		Cemented Stem (new)
		Porous Coated Stem (new)
	Trauma	Multifix LAG Screw
		Multifix Longnail
		Multifix Multi Screw
		Multifix Nail
Olympus Bio Materials	Total Knees	CW High Tibial Osteotomy - TriS Lateral HTO Blue (new)
		CW High Tibial Osteotomy - TriS Lateral HTO Gold (new)
		High Tibial Osteotomy - TriS (deleted)
		High Tibial Osteotomy - TriS Purple
		High Tibial Osteotomy - TriS Green
Omic Corporation	Trauma	OM Femoral Nail System 125°- Nail
		OM Femoral Nail System 130°- Nail
		OM Femoral Nail System - End Cap
		OM Femoral Nail System - Lag Screw
		OM Femoral Nail System - Locking Screw
		OM Femoral Nail System - Pin
Pega Medical	Trauma	FDTS-Distal Fixation with Pin and Femoral Proximal Fixation
		FDTS-Distal Fixation with Pin and Tibial Proximal Fixation
		FDTS-Distal Long Fixation and Femoral Proximal Fixation
		FDTS-Distal Short Fixation and Femoral Proximal Fixation
		FDTS-Distal Short Fixation and Tibial Proximal Fixation
Peter Brehm	Total Hips	ARR-Titan Acetabular Reconstruction Ring
		Ball head
		Bipolar head
		Vektor-Titan stems
Plus Orthopedics	Total Hips	SL-PLUS - Lateral
		SL-PLUS - Standard
Senko	Total Knees	MeraQ.uest. Knee System Femoral Component
		MeraQ.uest. Knee System Tibial Component Cemented
		MeraQ.uest. Knee System Tibial Insert+D1389
	Trauma	ALEXA NAIL SYSTEM 125°
		ALEXA NAIL SYSTEM 130°
Smith & Nephew	Total Hips	Anthology Hip System
		Basis 12/14 Taper
		Birmingham Hip Resurfacing – Resurfacing Head
		CPCS 12/14 Taper Stems
		Echelon Revision Hip System - Cemented Stem - 175mm
		Echelon Revision Hip System - Cemented Stem - 225mm
		Echelon Revision Hip System - Porous Bowed Stem
		Echelon Revision Hip System - Porous Straight Stem
		Echelon Stem - Porous Primary
		EPZ stem

Manufacturer	Category	Implant Overlay Template Family
		Image HA Hip System
		R3 Acetabular System
		REDAPT Revision Femoral Proximally Fluted Stems
		REDAPT Revision Femoral Sleeved Stems
		Reflection Cup - Cemented
		Reflection Cup - Porous/Coated
		Reflection Interfit Cup
		SL-PLUS Japan
		SL-PLUS MIA Hip Stem, Standard and Lateral
		SLR-PLUS Japan
		SMF Fixed Neck Surgical Templates
		SMF Modular Neck Implant Surgical Templates
		Spectron 12/14 Revision Taper Long Straight Stem
		Spectron 12/14 Revision Taper Neck Replacement Stem
		Spectrons Stem
		Synergy Select II Cemented Hip System
		Synergy Select II Porous Hip System
		Synergy Stem
		Synergy Stem – Cemented
		TANDEM Bipolar
	Total Knees	Bi-Cruciate Stabilized Knee System - Femoral
		Bi-Cruciate Stabilized Knee System - Tibial
		Genesis II Total Knee System – Femoral
		Genesis II Total Knee System Femoral Long Stem 100
		Genesis II Total Knee System Femoral Long Stem 150
		Genesis II Total Knee System Femoral Long Stem 200
		Genesis II Total Knee System Nonporous Tibial Baseplate
		Genesis II Total Knee System Porous Tibial
		Genesis II Total Knee System Tibial Alignment
		Genesis II Total Knee System Tibial Alignment Valgus
		Genesis II Total Knee System Tibial Long Stem 100
		Genesis II Total Knee System Tibial Long Stem 150
		Genesis II Total Knee System Tibial Long Stem 200
		Genesis II Total Knee System Tibial Wedge Angled Hemi
		Genesis II Total Knee System Tibial Wedge Full Plateau
		Genesis II Total Knee System Tibial Wedge Stepped Hemi
		JOURNEY II BCS OXINIUM FEMORAL (new)
		JOURNEY TIBIA BASEPLATE (new)
		JOURNEY II CR OXINIUM FEMORAL (new)
		Legion CR Cruciate Retaining Knee
		Legion CR Femoral Component Cemented (new)
		Legion Hinge Revision Knee System
		Legion Hinge Revision Knee System - Tibial
		Legion Oxinium CR Femoral Component (new)
		Legion Oxinium CR Narrow Femoral Component (new)
		Legion Oxinium PS Femoral Component (new)
		Legion Oxinium PS Narrow Femoral Component (new)
		Legion PS Femoral Component Cemented (new)
		Legion PS Posterior Stabilized Knee
		Legion Revision Knee System
		Legion Revision Knee System - Nonporous Tibial

Manufacturer	Category	Implant Overlay Template Family	
		Legion Revision Knee System - Offset Coupler	
		Legion Revision Knee System - Tibial Augment Full Angled	
		Legion Revision Knee System - Tibial Augment Full Plateau	
		Legion Revision Knee System - Tibial Augment Hemi Angled	
		Legion Revision Knee System - Tibial Augment Hemi Stepped	
		Profix Total Knee 4° Stem Solid	
		Profix Total Knee Conforming Plus Tibial Insert	
		Profix Total Knee Conforming Tibial Insert	
		Profix Total Knee Metaphyseal Tibial Stem	
		Profix Total Knee Neutral Stem Solid	
		Profix Total Knee Nonporous Tibial Base	
		Profix Total Knee Porous Femoral Component	
		Profix Total Knee Porous Tibial Base	
		Profix Total Knee System – Tibial Stem	
		Profix Total Knee System – Tibial Wedge	
		Solution EPP Femur	
		Solution EPP Revision Tibia	
		Solution EPP Tibia Basis	
		Total Shoulders	Cofield ² Total Shoulder – Eccentric Humeral Head
			Cofield ² Total Shoulder – Humeral Head
		Cofield ² Total Shoulder – Humeral Stem	
		Cofield ² Total Shoulder – Keeled All Poly Glenoid	
		Cofield ² Total Shoulder – Lateral Offset Humeral Head	
		Cofield ² Total Shoulder – Pegged All Poly Glenoid	
	Trauma	Trigen 1.5m Bow InterTAN Nail, 125° Ø10mm	
		Trigen 1.5m Bow InterTAN Nail, 125° Ø11.5mm	
		Trigen 1.5m Bow InterTAN Nail, 125° Ø13mm	
		Trigen 1.5m Bow InterTAN Nail, 130° Ø10mm	
		Trigen 1.5m Bow InterTAN Nail, 130° Ø11.5mm	
		Trigen 1.5m Bow InterTAN Nail, 130° Ø13mm	
		Trigen 2.0m Bow InterTAN Nail, 125° Ø10mm	
		Trigen 2.0m Bow InterTAN Nail, 125° Ø11.5mm	
		Trigen 2.0m Bow InterTAN Nail, 125° Ø13mm	
		Trigen 2.0m Bow InterTAN Nail, 130° Ø10mm	
		Trigen 2.0m Bow InterTAN Nail, 130° Ø11.5mm	
		Trigen 2.0m Bow InterTAN Nail, 130° Ø13mm	
		Trigen Adolescent Trochanteric Antegrade Nail, 130° Ø8.5mm	
		Trigen Antegrade Tibial Only, Ø8.5mm	
		Trigen Femoral Antegrade Nail, Ø10mm	
		Trigen Femoral Antegrade Nail, Ø11.5mm	
		Trigen Femoral Antegrade Nail, Ø13mm	
		Trigen Femoral Antegrade Nail, Ø14.5mm	
	Trigen Femoral Antegrade Nail, Ø16mm		
	Trigen Long Bent Humeral Nail		
	Trigen Meta Retrograde Femoral Nail, Ø10mm		
	Trigen Meta Retrograde Femoral Nail, Ø11.5mm		
	Trigen Meta Retrograde Femoral Nail, Ø13mm		
	Trigen Meta Tibia Nail, Ø10mm		
	Trigen Meta Tibia Nail, Ø11.5mm		
	Trigen Meta Tibia Nail, Ø13mm		
	Trigen Meta Tibia Nail, Ø8.5mm		

Manufacturer	Category	Implant Overlay Template Family
		Trigen Proximal Bent Humeral Nail
		Trigen Proximal Straight Humeral Nail
		Trigen Retrograde Femoral or Antegrade Tibia, Ø10mm
		Trigen Retrograde Femoral or Antegrade Tibia, Ø11.5mm
		Trigen Retrograde Femoral or Antegrade Tibia, Ø13mm
		Trigen Retrograde Supracondylar Nail, Ø11mm
		Trigen Retrograde Supracondylar Nail, Ø13mm
		Trigen Short Trochanteric Antegrade Nail
		Trigen Trochanteric Antegrade Nail, 130° Ø10mm
		Trigen Trochanteric Antegrade Nail, 130° Ø11.5mm
		Trigen Trochanteric Antegrade Nail, 130° Ø13mm
		Trigen Trochanteric Antegrade Nail, 135° Ø10mm
		Trigen Trochanteric Antegrade Nail, 135° Ø11mm
		Trigen Trochanteric Antegrade Nail, 135° Ø13mm
		Stryker
Solar Total Elbow System – Ulnar Component		
Total Hips	ABGII Acetabular System Ceramic Insert	
ABGII Acetabular System Hooded PE Insert		
ABGII Acetabular System Standard PE Insert		
ABGII Hip Stem		
Accolade C Cemented Hip Stem		
Accolade Hfx Hip Stem		
Accolade II Femoral Hip Stem		
ACCOLADE TMZF Hip Stem		
Anato Hip System Anteverted Stem		
Anato Hip System Neutral Stem		
Austin Moore Hip Stem Fenestrated Narrow - Curved Stem		
Austin Moore Hip Stem Fenestrated Regular - Curved Stem		
CentPillar HA Stem		
CentPillar TMZF Stem - Primary TMZF/Plasma Spray HA Hip Stem		
Centrax Bipolar System Centrax Bipolar Head		
Citation TMZF Hip Stem		
Crossfire All-Poly Cup		
Definition PM Hip Stem		
Exeter Contemporary Flanged Cup		
Exeter Hooded Contemporary Acetabular Cup		
Exeter Total Hip High Profile Acetabular Cup		
Exeter Total Hip System		
Exeter Total Hip System Long Stems		
Exeter Total Hip X-tra Cover High Profile Acetabular Cup		
Exeter Total Hip X-tra Cover Low Profile Acetabular Cup		
EXETER X3 RimFit cup		
HIPSTAR Hip Stem		
Meridian TMZF Hip Stem		
MITCH TRH System Acetabular Cup		
MITCH TRH System Proximal Epiphyseal Replacement		
MITCH TRH System Resurfacing Head		
MITCH TRH System V40 Modular Head - Offset +4		
MITCH TRH System V40 Modular Head - Offset +8		
MITCH TRH System V40 Modular Head - Offset 0		
MITCH TRH System V40 Modular Head - Offset -4		

Manufacturer	Category	Implant Overlay Template Family
		Mobile Bearing Hip System - Modular Dual Modality Liner
		Omnifit EON Cemented Hip Stem
		Omnifit HFX Cemented Hip Stem
		Rejuvenate SPT Modular Stem
		Reliance PF Hip Stem
		Restoration ADM Acetabular Shell
		Restoration C Cemented Revision Hip Stem
		Restoration C Graft Tube
		Restoration HA Revision Hip System – 205 mm Calcar Replacement Stem
		Restoration HA Revision Hip System – Bowed Stem
		Restoration HA Revision Hip System – Straight Stem
		Restoration Modular Broached Body
		Restoration Modular Calcar Body
		Restoration Modular Cone Body
		Restoration Modular Conical Distal Stem (Bowed)
		Restoration Modular Conical Distal Stem (Straight)
		Restoration Modular Fluted Distal Stem
		Restoration Modular Milled Body
		Restoration Modular MT3 Body
		Restoration Modular Plasma Distal Stem (Bowed)
		Restoration Modular Plasma Distal Stem (Straight)
		Restoration PS Revision Hip System – Bowed Fully Coated Long Stem
		Restoration PS Revision Hip System – Bowed Long Stem
		Restoration PS Revision Hip System – Bowed Stem
		Restoration PS Revision Hip System – Straight Stem
		Secur-Fit Advanced Femoral Hip Stem, 127 deg
		Secur-Fit Advanced Femoral Hip Stem, 132 deg
		Secur-Fit HA Hip Stem
		Secur-Fit HA Max Hip Stem
		Secur-Fit Plus Hip Stem
		Secur-Fit Plus Max Hip Stem
		Super EON Stem – Primary Hip Stem, 127deg
		Super EON Stem – Primary Hip Stem, 132deg
		Super Secur-Fit – Primary HA Coated Hip Stem, 127deg
		Super Secur-Fit – Primary HA Coated Hip Stem, 132deg
		Super Secur-Fit Plus HA Stem
		Symax Cementless Hip Stem
		TriAD HA Hemispherical Cup
		TriAD PSL HA Cluster Cup
		TriAD PSL HA Cup, Solid-Back
		Trident Acetabular System - All Poly Constrained Acetabular Insert
		Trident Alumina Insert
		Trident Crossfire Polyethylene Insert
		Trident ECCENTRIC Crossfire Polyethylene Insert
		Trident ECCENTRIC X3 Polyethylene Insert
		Trident ELEVATED RIM Crossfire Polyethylene Insert
		Trident ELEVATED RIM X3 Polyethylene Insert
		Trident Hemispherical Cup
		Trident PSL Cup
		Trident X3 Polyethylene Insert
		Tritanium Primary Hemispherical Shell, Cluster Hole

Manufacturer	Category	Implant Overlay Template Family
		Tritanium Primary Hemispherical Shell, Solidback
	Total Knees	AVON Patello-Femoral Arthroplasty
		Duracon Total Knee System – Universal Baseplates
		Duracon Total Knee System Standard Femoral Component
		Duracon Total Knee System Condylar Tibial Insert
		Duracon Total Knee System Cruciform Tibial Baseplate Beaded
		Duracon Total Knee System Full Posterior Spacers
		Duracon Total Knee System Half Posterior Spacers
		Duracon Total Knee System Press Fit Stem Extenders
		Duracon Total Knee System Tibial Half Augment
		Duracon Total Knee System Total Stabilizer Femoral Components
		Duracon Total Knee System Total Stabilizer Tibial Inserts
		Duracon Total Knee System TS Femoral Augmentation Distal Blocks
		GMRS Knee System All Poly Tibial Component
		GMRS Knee System Cemented Stem (Bowed with Body)
		GMRS Knee System Cemented Stem (Bowed with Body) 6.58"
		GMRS Knee System Cemented Stem (Bowed with Body) 7.58"
		GMRS Knee System Cemented Stem (Bowed without Body)
		GMRS Knee System Cemented Stem (Bowed without Body) 5.44"
		GMRS Knee System Cemented Stem (Bowed without Body) 6.44"
		GMRS Knee System Cemented Stem (Straight with Body)
		GMRS Knee System Cemented Stem (Straight)
		GMRS Knee System Connection Pieces
		GMRS Knee System Extension Pieces
		GMRS Knee System Femoral Component
		GMRS Knee System Press-Fit Stem (Bowed)
		GMRS Knee System Press-Fit Stem (Straight Fluted)
		GMRS Knee System Proximal Tibia Component
		GMRS Knee System Tibial Inserts
		GMRS Proximal Femoral Component
		Kinemax Total Knee System Condylar LS Insert
		Kinemax Total Knee System Femoral Component - Condylar
		Kinemax Total Knee System Modular Tibial Baseplate Non-Beaded
		MRH Knee System Femoral Component
		MRH Knee System Tibial Baseplate
		MRH Knee System Tibial Insert
		Scorpio CR Femoral Component - Basic
		Scorpio CR Tibial Insert
		Scorpio NRG Total Knee System - CR Femoral Component
		Scorpio NRG Total Knee System - CR Inserts
		Scorpio NRG Total Knee System - HA Coated CR Femoral Component
		Scorpio NRG Total Knee System - HA Coated PS Femoral Component
		Scorpio NRG Total Knee System - PS Femoral Component
		Scorpio NRG Total Knee System - PS Inserts
		Scorpio SuperFlex – CR HA Femoral Component
		Scorpio SuperFlex – CR Inserts
		Scorpio SuperFlex – CR Micro Structure Femoral Component
		Scorpio SuperFlex – CR Waffle Femoral Component
		Scorpio SuperFlex – PS HA Femoral Component
		Scorpio SuperFlex – PS Inserts
		Scorpio SuperFlex – PS Micro Structure Femoral Component

Manufacturer	Category	Implant Overlay Template Family
		Scorpio SuperFlex – PS Waffle Femoral Component
		Scorpio Total Knee System – TS Constrained Tibial Bearing Insert
		Scorpio Total Knee System Cemented Stem Extenders
		Scorpio Total Knee System Fluted Stem Extenders
		Scorpio Total Knee System TS Femoral Augmentation Posterior Block
		Scorpio Total Knee System TS Femoral Component
		Scorpio Total Knee System TS Offset Adapter
		Scorpio Total Knee System TS Tibial Augmentation Full 10 mm Block
		Scorpio Total Knee System TS Tibial Augmentation Full 5 deg Wedge
		Scorpio Total Knee System TS Tibial Augmentation Half Block
		Scorpio Total Knee System TS Tibial Tray
		Series 7000 Primary Tibial Tray
		Series 7000 Total Knee – Primary Tibial Tray HA
		Series 7000 Total Knee – Primary Tibial Tray Microstructured
		Series 7000 Total Knee – Primary Tibial Tray Waffled
		Triathlon Knee System – PS Tibial Insert
		Triathlon Knee System – TS Tibial Insert
		Triathlon Knee System Primary Tibial Tray
		Triathlon Knee System PS Femoral Component
		Triathlon Knee System TS Cemented Stem
		Triathlon Knee System TS Femoral Component
		Triathlon Knee System TS Offset Adapter
		Triathlon Knee System TS Stem Extender
		Triathlon Knee System TS Tri Fluted Cementless Stem
		Triathlon Total Knee System CR Femoral Component Cemented w/Lugs
		Triathlon Total Knee System CR Femoral Component Cementless w/Lugs
		Triathlon Total Knee System CR Tibial Insert
		Triathlon Total Knee System TS Femoral Distal Augment
		Triathlon Total Knee System TS Femoral Posterior Augment
		Triathlon TS Tibial Baseplate
	Total Shoulders	ReUnion TSA System - Glenoid Component - Keeled
		ReUnion TSA System - Glenoid Component - Pegged
		ReUnion TSA System - Humeral Stem Component - Cemented
		ReUnion TSA System - Humeral Stem Component - Cemented Long
		ReUnion TSA System - Humeral Stem Component - Press-fit
		Solar Total Shoulder Bipolar Heads
		Solar Total Shoulder Bipolar Modular Heads
		Solar Total Shoulder Dual Offset Humeral Head
		Solar Total Shoulder Humeral Component
		Solar Total Shoulder Humeral Head
		Solar Total Shoulder System - Glenoid
	Trauma	AxSOS Locking Screw
		AxSOS Locking Screw, Sterile packaged product
		AxSOS Periprosthetic Locking Screw
		AxSOS Periprosthetic Locking Screw, Sterile packaged product
		Gamma3 End Cap Std, Ti
		Gamma3 Lag Screw, Ti
		Gamma3 Long Nail Kit R1.5, Ti, Ø 10 120°
		Gamma3 Long Nail Kit R1.5, Ti, Ø 10 125°
		Gamma3 Long Nail Kit R1.5, Ti, Ø 10 130°
		Gamma3 Long Nail Kit R2.0, Ti, Ø 11 120°

Manufacturer	Category	Implant Overlay Template Family
		Gamma3 Long Nail Kit R2.0, Ti, Ø 11 125°
		Gamma3 Long Nail Kit R2.0, Ti, Ø 11 130°
		Gamma3 Long Nail Kit R2.0, Ti, Ø 13 120°
		Gamma3 Long Nail Kit R2.0, Ti, Ø 13 125°
		Gamma3 Long Nail Kit R2.0, Ti, Ø 13 130°
		Gamma3 Trochanteric Nail Kit, Ti-120
		Gamma3 Trochanteric Nail Kit, Ti-125
		Gamma3 Trochanteric Nail Kit, Ti-130
		Proximal Humerus Plate
		Proximal Humerus Plate, Sterile packaged product
		T2 Advanced Compression Screw, Femur
		T2 Compression Screw for Tibial Nail
		T2 Condyle Screw Nut, Ti
		T2 Condyle Screw Ø 5
		T2 End-Cap for Femur and Tibia
		T2 End-Cap for Femur, Tibia and Reconstruction Nails
		T2 End-Cap for Humerus
		T2 End-Cap for Proximal Humeral Nail
		T2 End-Cap for Reconstruction Nail
		T2 Femoral Nail A/R 1500R
		T2 Femoral Nail A/R Ø10mm
		T2 Femoral Nail A/R Ø11mm
		T2 Femoral Nail A/R Ø12mm
		T2 Femoral Nail A/R Ø13mm
		T2 Femoral Nail A/R Ø14mm
		T2 Femoral Nail A/R Ø15mm
		T2 Femoral Nail A/R Ø8mm
		T2 Femoral Nail A/R Ø9mm
		T2 Humeral Nail
		T2 Knee Arthrodesis Nail Ø10mm
		T2 Knee Arthrodesis Nail Ø11.5mm
		T2 Knee Arthrodesis Nail Ø13mm
		T2 Lag Screw for Reconstruction Nail
		T2 Locking Screw, Fully Threaded Ø4mm
		T2 Locking Screw, Fully Threaded Ø5mm
		T2 Locking Screw, Shaft Screw
		T2 Proximal Humeral Nail, Long
		T2 Proximal Humeral Nail, Short
		T2 Reconstruction Nail 1500R_size10mm
		T2 Reconstruction Nail 1500R_size11mm
		T2 Reconstruction Nail 1500R_size12mm
		T2 Reconstruction Nail 1500R_size13mm
		T2 Reconstruction Nail 1500R_size9mm
		T2 Set Screw for Reconstruction Nail, Antegrade
		T2 Set Screw for Reconstruction Nail, Recon
		T2 Supracondylar Nail Ø10mm
		T2 Supracondylar Nail, Long
		T2 Supracondylar Nail, Short
		T2 Tibial Nail , Distal Tibia
		T2 Tibial Nail, Standard_10mm
		T2 Tibial Nail, Standard_11mm

Manufacturer	Category	Implant Overlay Template Family
		T2 Tibial Nail, Standard_12mm
		T2 Tibial Nail, Standard_13mm
		T2 Tibial Nail, Standard_8mm
		T2 Tibial Nail, Standard_9mm
		T2 Washer Round
		T2 Washer Square
	Unicondylar Knees	ELIUS Unicompartamental Knee System Femoral Component Implants
		Triathlon Partial Knee Resurfacing (PKR) Femoral Component - Cemented
		Triathlon Partial Knee Resurfacing (PKR) Tibial Baseplate - Cemented
Synthes	Trauma	1.3mm Cortex Screw, Self-Tapping - Stainless Steel
		1.3mm Cortex Screw, Self-Tapping - Titanium
		1.3mm Extended H-Plate
		1.3mm Straight Plate
		1.3mm T-Plate
		1.3mm Y-Plate
		1.5mm Cortex Screw
		1.5mm Extended H-Plate
		1.5mm Straight Plate
		1.5mm T-Plate 3 Hole Head
		1.5mm T-Plate 4 Hole Head
		1.5mm Y-Plate
		1.8mm Buttress Pin
		10mm TI Cannulated Distal Femoral Nail
		10mm TI Cannulated Retrograde Femoral Nail-EX
		11mm TI Cannulated Retrograde Femoral Nail-EX
		12mm TI Cannulated Distal Femoral Nail
		12mm TI Cannulated Retrograde Femoral Nail-EX
		13mm TI Cannulated Distal Femoral Nail
		13mm TI Cannulated Retrograde Femoral Nail-EX
		14mm TI Cannulated Retrograde / Antegrade Femoral Nail-EX
		14mm TI Cannulated Retrograde Femoral Nail-EX
		14mm TI Expert A2FN Cannulated
		15mm TI Cannulated Retrograde / Antegrade Femoral Nail-EX
		15mm TI Cannulated Retrograde Femoral Nail-EX
		2.0mm Cortex Screw
		2.0mm Straight Plate
		2.0mm T-Plate
		2.0mm Y-Plate
		2.4/2.7mm Locking X-Plate ,Stainless Steel
		2.4/2.7mm Locking X-Plate ,Titanium
		2.4mm Cortex Screw Self-Tapping
		2.4mm Cortex Screw Self-Tapping with T8 Stardrive Recess
		2.4mm LCP Distal Radius L Plate -90°2H HD
		2.4mm Locking Screw Self-Tapping with StarDrive Recess
		2.4mm Straight Plate
		2.4mm TI LCP Distal Radius Plate - Straight
		2.4mm TI LCP Distal Radius Plate - Straight / Long
		2.4mm T-Plate
		2.4mm Y-Plate
		2.4mmVA-LCP Two Column Distal Radius Volar Plate ,Stainless Steel
		2.4mmVA-LCP Two Column Distal Radius Volar Plate ,Titanium

Manufacturer	Category	Implant Overlay Template Family
		2.4mmVA-LCP Two Column Distal Radius Volar Plate, Narrow ,Stainless Steel
		2.4mmVA-LCP Two Column Distal Radius Volar Plate, Narrow ,Titanium
		2.7mm Cortex Screw - Hex Drive
		2.7mm Locking Screw Self-Tapping with StarDrive Recess
		3.0mm Cannulated Screw, Short Thread
		3.5mm Cortex Screw - Hex Drive
		3.5mm Curved Reconstruction Plate
		3.5mm LCP Medial Distal Humerus Plate
		3.5mm LCP Metaphyseal Plate
		3.5mm LCP Olecranon Plate - Left
		3.5mm LCP Olecranon Plate - Right
		3.5mm LCP Plate
		3.5mm LCP Posterolateral Distal Humerus Plate
		3.5mm LCP Posterolateral Distal Humerus Plate - Lateral Support
		3.5mm LCP Proximal Tibia Plate, Right
		3.5mm LCP T-Plate 3 hole Head Oblique - Right
		3.5mm LCP T-Plate 3 Hole Head Right Angle
		3.5mm LCP T-Plate 4 Hole Head - Right Angle
		3.5mm Locking Screw Self Tapping Stardrive
		3.5mm Pediatric LCP Hip Plate
		3.5mm Reconstruction Plate
		3.5mm Ti LCP Metaphyseal Plate
		3.5mm T-Plate, Oblique, 3 Hole Head
		3.5mm T-Plate, Right Angle
		3.9mm Titanium Locking Bolt
		4.0mm Cancellous Bone Screw, Fully Threaded
		4.0mm Cancellous Bone Screw, Partially Threaded
		4.0mm Cannulated Screw, Short Thread, Self-Drilling, Self-Tapping, Stainless Steel
		4.0mm Cannulated Screw, Short Thread, Self-Drilling, Self-Tapping, Titanium
		4.0mm Ti Locking Screw Stardrive For IM Nails
		4.5mm Cannulated Screw, Partially Threaded
		4.5mm Cortex Screw with Large Hex
		4.5mm Cortex Screw, Self Tapping
		4.5mm Curved Broad LCP Plate
		4.5mm LCP T-Plate
		4.5mm Malleolar Screw
		4.5mm Nut
		4.5mm Shaft Screw
		4.9mm Titanium Locking Bolt
		5.0mm Locking Screw, Self Tapping
		5.0mm Pediatric LCP Hip Plate
		5.0mm Periprosthetic Locking Screw, Self-Tapping, Stainless Steel
		5.0mm Periprosthetic Locking Screw, Self-Tapping, Titanium
		5.0mm Ti Locking Screw Stardrive For IM Nails
		6.0mm TI Locking Screw w-T25 Stardrive for IM Nails
		6.5mm Cancellous Bone Screw - 16mm Thread
		6.5mm Cancellous Bone Screw - 32mm Thread
		6.5mm Cancellous Bone Screw Fully Threaded
		6.5mm Cannulated Screw, 16mm Thread, Self-Drilling, Self-Tapping, Stainless Steel
		6.5mm Cannulated Screw, 16mm Thread, Self-Drilling, Self-Tapping, Titanium
		6.5mm Cannulated Screw, 32mm Thread, Self-Drilling, Self-Tapping, Stainless Steel

Manufacturer	Category	Implant Overlay Template Family
		6.5mm Cannulated Screw, 32mm Thread, Self-Drilling, Self-Tapping, Titanium
		7.3mm Cannulated Screw 16mm Thread
		7.5mm TI Solid Humeral Nail
		90 Degree Child Osteotomy Plate
		90 Degree Infant Osteotomy Plate
		9mm TI Cannulated Distal Femoral Nail
		9mm TI Cannulated Proximal Humeral Nail-EX
		9mm TI Cannulated Retrograde Femoral Nail-EX
		Cerclage Positioning Pin Stainless Steel
		Cerclage Positioning Pin Titanium
		DCS Plate - 95 degrees
		DHS/DCS - 130 degree DHS Plate - Short Barrel
		DHS/DCS - 130 degree DHS Plate - Standard Barrel
		DHS/DCS - 135 degree DHS Plate - Standard Barrel
		DHS/DCS - 140 degree DHS Plate - Short Barrel
		DHS/DCS - 140 degree DHS Plate - Standard Barrel
		DHS/DCS - 145 degree DHS Plate - Short Barrel
		DHS/DCS - 145 degree DHS Plate - Standard Barrel
		DHS/DCS - 150 degree DHS Plate - Short Barrel
		DHS/DCS - 150 degree DHS Plate - Standard Barrel
		DHS/DCS Compression Screw
		DHS/DCS Lag Screw
		End Cap For Ti Distal Femoral Nails
		Expert Tibial Nail with Proximal Bend
		Large Fragment - 4.5mm Broad LCP Plate
		Large Fragment - 4.5mm LCP Plate
		Large Fragment - 4.5mm LCP Straight Reconstruction Plate
		Large Fragment - 4.5mm Ti Narrow LCP Plate
		LCP Anterolateral Distal Tibial Plate 3.5 - Stainless Steel
		LCP Anterolateral Distal Tibial Plate 3.5 - Titanium
		LCP Distal Femur Plates (updated)
		LCP Distal Radius Plate Volar - Left
		LCP Distal Tibial Plate 2.7/3.5, medial
		LCP Lateral Fibula Plates ,Stainless Steel
		LCP Lateral Fibula Plates ,Titanium
		LCP Posterolateral Fibula Plates ,Stainless Steel
		LCP Posterolateral Fibula Plates ,Titanium
		LCP Proximal Humerus Plates - Long
		LCP Proximal Humerus Plates - Standard
		LCP Proximal Lateral Tibia Plate
		LCP Reconstruction Plates
		LCP Volar Distal Radius Plate Extra - Articular / Left
		LCP Volar distal Radius Plate Extra - Articular / Right
		LCP Volar Distal Radius Plate Extra-Articular 4H Head-Left
		LCP Volar Distal Radius Plate Extra-Articular 4H Head-Right
		One-Third Tubular Plate, with Collar
		PFNA Proximal Femoral Nails
		PFNA-II Blades Enhanced
		PFNA-II Proximal Femoral Nails
		Recon Screw-LFN
		Small Fragment - 3.5mm LCP Medial Distal Tibia Plate - Left

Manufacturer	Category	Implant Overlay Template Family		
		Small Fragment - 3.5mm LCP Medial Distal Tibia Plate - Right		
		Small Fragment - 3.5mm LCP Proximal Tibia Plate-Left		
		Small Fragment - 3.5mm Ti LCP Plate		
		Small Fragment - LCP Distal L Plate +90 degree 2H HD		
		Small Fragment - LCP Distal L Plate -90 degree 2H HD		
		Small Fragment - LCP Proximal Humerus Plate - Long		
		Small Fragment - LCP Proximal Humerus Plate - Standard		
		Small Fragment - LCP Reconstruction Plate		
		Spherical Washers, for 6.5mm & 7.3mm Cannulated Screw		
		Threaded Guide Wire		
		Ti Cannulated Humeral Nail-Ex		
		Ti Cannulated Retrograde Antegrade Femoral Nail-Ex		
		Ti Distal Femur LISS Plate		
		TI Dorsal Distal Radius Plate		
		Ti End Cap T40 Stardrive Ext for Femoral Nails-EX		
		Ti Expert A2FN Cannulated		
		Ti Expert A2FN End Cap Cannulated Extension		
		Ti LCP Distal Radius Plate Volar - Right		
		Ti Proximal Tibia LISS Plate		
		Ti Spiral Blade for Humeral Nail		
		TI Spiral Blade For TI Retrograde Femoral Nails-EX		
		Ti Tomofix Lateral Distal Femur Plate		
		Ti Tomofix Lateral High Tibia Plate		
		Ti Tomofix Medical High Tibia Plate		
		Ti Volar Distal Radius Plate		
		Titanium End Cap		
		Titanium End Cap For Titanium Solid Tibial Nails		
		TomoFix Medial Distal Femur Version 2011		
		TomoFix Medial High Tibia		
		VA Cannulated Locking Screw ø5.0mm		
		VA Locking 2.4/2.7mm Screw Stardrive, Self-Tapping, Stainless Steel		
		VA Locking 2.4/2.7mm Screw Stardrive, Self-Tapping, Titanium		
		VA Locking Screw Stardrive ø5.0mm, self-tapping		
		VA Periprosthetic Locking Screw Stardrive ø5.0mm, self-tapping		
		VA-LCP Condylar Plate 4.5/5.0, Stainless Steel		
		VA-LCP Condylar Plate 4.5/5.0, Titanium		
		Washers		
		Teijin Nakashima Medical	Ankle	FINE Total Ankle System Talus Component
				FINE Total Ankle System Tibial Socket./Tibial Socket Holder
			Total Elbows	K-NOW Total Elbow System (Humeral Component/Humeral Stem)
				K-NOW Total Elbow System (Ulnar Component/Ulnar Stem)
K-NOW+C Total Elbow System (Humeral Component/Humeral Stem)				
K-NOW+C Total Elbow System (Ulnar Component/Ulnar Stem)				
K-NOW+C Total Elbow System (Ulnar Component-All Polyethylene Type)				
Total Hips	4-U Stem			
	BP Cup(22)			
	BP Cup(26)			
	Delta-Lock Stem 170mm			
	Delta-Lock Stem 225mm			
	Delta-Lock Stem 280mm			
	Dolphin CLF cup			

Manufacturer	Category	Implant Overlay Template Family	
		Dolphin CLH cup	
		Dolphin Stem	
		F.S. Stem	
		Joinus Hip Stem	
		Teijin Nakashima Medical THA Cup 3 hole (Elevate Type)	
		Teijin Nakashima Medical THA Cup 3 hole (Flat Type)	
		Teijin Nakashima Medical THA Cup Multi hole (Elevate Type)	
		Teijin Nakashima Medical THA Cup Multi hole (Flat Type)	
		Preserve Alpha (updated)	
		Preserve Beta (updated)	
		Prime-N Collared 130°	
		Prime-N Collared 135°	
		Prime-N Collarless 130°	
		Prime-N Collarless 135°	
		TNH Cup Elevate	
		TNH Cup Flat	
		Twinkle cup	
		Total Knees	Fine Total Knee System CR Femoral Component
			Fine Total Knee System CR Tibial Component
			Fine Total Knee System CR-Cementless Femoral Component
		Fine Total Knee System CR-Cementless Tibial Component	
		Fine Total Knee System PS Femoral Component	
		Fine Total Knee System PS Tibial Component	
		FNK Total Knee System CR Femoral Component	
		FNK Total Knee System CR Tibial Component	
		FNK Total Knee System PS Femoral Component	
		FNK Total Knee System PS Tibial Component	
	Trauma	N-Radial Head Plate	
		ONI ELBOW SYSTEM	
	Tornier	Ankle	Salto Talaris Fixed Bearing
			Salto Talaris Talar Implant
			Salto Talaris Tibial Implant
		Total Hips	Cemented Dual Mobility Acetabular Cup
Cementless Dual Mobility Acetabular Cup			
Cementless Dual Mobility Acetabular Cups With Flanges And Hook			
Cementless Dual Mobility Cup with flanges			
Total Shoulders		Aequalis ASCEND Humeral Head	
		Aequalis ASCEND Humeral Stem (Non Coated)	
		Aequalis ASCEND PTC Humeral Stem (Coated)	
		Aequalis Head CoCr	
		Aequalis PerFORM - CortiLoc Glenoid	
		Aequalis PerFORM - Keeled Glenoid	
		Aequalis PerFORM - Pegged Glenoid	
		Aequalis Standard Stem	
		Humeral Head - Inspyre	
UNITECH (new)		Arm (new)	Stellar2 (new)
Wright Medical	Foot	SWANSON flexible hinge toe	
	Hand	SWANSON flexible hinge finger joint	
	Total Hips	ANCA Fit Cup Polyethylene Insert 10° Offset	
		ANCA Fit Cup Standard Insert	
ANCA Fit Hip Stem			

Manufacturer	Category	Implant Overlay Template Family		
		Bipolar Cup		
		Conserve Cup		
		Conserve Femoral Component		
		Conserve Plus / Total 10mm Acetabular Shells		
		Conserve Plus / Total 6mm Acetabular Shells		
		Conserve Plus Femoral Component		
		Conserve Plus SF		
		Conserve Total Super-Fix Acetabular Cup		
		Dynasty® Acetabular Shell		
		Lineage Acetabular Cup Deep Profile		
		Lineage Acetabular Cup Standard Profile		
		Perfecta Plasma Spray II		
		Perfecta RS II SLT Taper		
		Perfecta SLT Taper Short Neck		
		Perfecta SLT Taper Slim Neck		
		Profemur Hip System Distal Stem Cylindrical Plasma Spray		
		Profemur Hip System Distal Stem Cylindrical Spline with Slot		
		Profemur Hip System Distal Stem Tapered Splined		
		Profemur Hip System Modular Neck		
		Profemur Hip System Proximal Body		
		Profemur Renaissance Hip Stem		
		Profemur Tapered Hip Stem		
		Profemur Xm		
		ProfemurE Hip Stem		
		ProfemurL Hip Stem		
		ProfemurTL Hip Stem		
		ProfemurZ Hip Stem		
		Total Knees		Advance Offset Tibial Base
				Advance Revision Knee System Femoral Component
				Advance Revision Knee System Tibial Component
				Advance Stemmed Medial Pivot Knee System
				Advance Total Knee System All-Poly Universal Tibial Component Posterior Stabilized
	Advance Total Knee System Fixed Stem Tibial Base Posterior Stabilized Insert			
	Advance Total Knee System Posterior Stabilized Universal Femoral Component			
	Advance Total Knee System Primary Femoral Implant			
	Advance Total Knee System Tibial Base And Insert			
	EVOLUTION Medial Pivot Total Knee System CS/CR Primary Femoral Component			
	EVOLUTION Medial Pivot Total Knee System Primary Cemented Base with CS Insert or CR Insert			
	Trauma		MICRONAIL	
	Zimmer	Ankle	TM Total Ankle - Prolong HXPE Insert	
			TM Total Ankle - Talar Component	
			TM Total Ankle - Tibial Base	
		Hand	Elogenics Finger - Basis	
			Elogenics Finger - Caput	
		Total Elbows		C/M Humeral Component
				C/M Humeral Component Long Flange
				C/M Ulna Component
				C/M Ulna Component Long Flange

Manufacturer	Category	Implant Overlay Template Family
		Nexel Humeral Prosthesis
		Nexel Ulnar Prosthesis
	Total Hips	A2 Primary Stem
		ACA Stem
		Alloclassic Variall Cup
		Alloclassic Variall SLV Stem
		Alloclassic Zweymueller CSF Shell
		Alloclassic Zweymueller SL Stem
		Alloclassic Zweymueller SL-64 Stem
		Alloclassic Zweymueller SLA Stem
		Alloclassic Zweymueller SL-HAC Stem
		Alloclassic Zweymueller SLL Revision Stem
		Alloclassic Zweymueller SLL Revision Stem
		Alloclassic Zweymueller SLO Stem
		Allofit Alloclassic Shell
		Allofit IT S Shell
		Allofit IT Shell
		Allofit-S Alloclassic Shell
		Alpha Cerasul Insert neutral liner
		Alpha Durasul Insert hooded liner
		Alpha Metasul Insert hooded liner
		Alpha Metasul Insert neutral
		Alumina Femoral Head 12/14
		APR Acetabular Shell, Protrusio
		APR Acetabular Shell, Standard
		APR CoCr Nonporous Femoral
		APR Collarless Femoral Large Body Distally Textured
		APR Collarless Femoral Standard Body Distally Textured
		APR Femoral Large Body Distally Textured
		APR Femoral Large Body Fully Textured
		APR Femoral Oversized Distally Textured
		APR Femoral Oversized Fully Textured
		APR Femoral Standard Body Distally Textured
		APR HA Femoral Large Body Distally Textured
		APR HA Femoral Oversized Distally Textured
		APR HA Femoral Standard Body Distally Textured
		APS Natural-Hip Stem Anteverted
		APS Natural-Hip Stem RM Non-Anteverted
		APS Natural-Hip Stem Standard Non-Anteverted
		Avenir Hip Stem Cemented - Lateral
		Avenir Hip Stem Cemented - Standard
		Avenir-Mueller Stem Lateral
		Avenir-Mueller Stem Standard
		BEO Stem
		BIOLOX delta Ceramic Femoral Head 12/14
		BIOLOX delta Taper Liner
		BIOLOX forte Ceramic Femoral Head 12/14
		BIOLOX OPTION Femoral Head / Adapter 12/14
		BIOLOX OPTION Femoral Head / Adapter 8/10
		Bipolar Assembly CoCr
Bipolar Metal Shell – with Poly Liner Assembly – 22, 26, 28 & 32mm		

Manufacturer	Category	Implant Overlay Template Family
		Bipolar Metal Shell – with Poly Liner Assembly – 22mm
		Bipolar Metal Shell – with Poly Liner Assembly – 26, 28 & 32mm
		Bipolar Metal Shell – with Poly Liner Assembly – 26mm
		Bipolar Metal Shell – with Poly Liner Assembly – 28mm
		Bipolar Metal Shell – with Poly Liner Assembly – 32mm
		Brunswick Cup
		Burch-Schneider Cage
		Burch-Schneider Reinforcement Cage
		Cerasul Ceramic Femoral Head 12/14
		Cerasul Ceramic Femoral Head 8/10
		CLS Brevius Stem with Kinectiv Technology
		CLS Spotorno 125° Stem
		CLS Spotorno 135° Stem
		CLS Spotorno 145° Stem
		CLS Spotorno 145° Stem with HAC
		CLS Spotorno Expansion Shell
		CoCr Femoral Head 12/14
		CoCr Slim
		Cone Prosthesis
		Continuum Shell, Cluster-Hole
		Continuum Shell, Multi-Hole
		Continuum Shell, Uni-Hole
		Converge CSTi Cluster-Hole Porous Shell with Sealed Screw Holes
		Converge CSTi Hemispherical Porous Shell with Sealed Screw Holes
		Converge CSTi Porous Multi Hole Shell
		Converge CSTi Rim Flare Porous Shell
		Converge Poly Liner Durasul
		Converge Protrusio Porous Shell with Sealable Screwholes
		Converge Rim Flare Porous Shell with Spikes and Sealed Screw Holes
		CPT Primary Ext Offset
		CPT Primary Extra Ext Offset
		CPT Primary Small
		CPT Primary Std Offset
		CPT Revision Long
		Durasul Bipolar CoCr Shell
		Durasul CLS Spotorno Insert
		Durasul CoCr Femoral Head 12/14
		Durasul CoCr Femoral Head 8/10
		Durasul PE Alpha Insert
		Durom Hip System - Acetabular Component
		Durom Hip System - Femoral Component
		Durom Hip System - US Acetabular Component
		Emeraude Stem cemented
		Emeraude Stem uncemented
		Epoch Stem
		Exafit Dysplasia Stem
		Exafit Revision Stem
		Exafit Stem
		Femoral Head 6°
		Fitmore Hip Stem A
		Fitmore Hip Stem B

Manufacturer	Category	Implant Overlay Template Family
		Fitmore Hip Stem B Ext. Offset
		Fitmore Hip Stem C
		Fitmore Shell, with fins
		Fitmore Shells, with Screw Cones
		HGP II Porous Acetabular Cup
		JB Stem
		Lamella Cup
		Lamella Stem
		Longevity Highly Crosslinked IT Liner, Elevated
		Longevity Highly Crosslinked IT Liner, Neutral
		Longevity Highly Crosslinked Polyethylene Liner, 10° Elevated
		Longevity Highly Crosslinked Polyethylene Liner, 20° Elevated
		Longevity Highly Crosslinked Polyethylene Liner, 7mm Offset
		Longevity Highly Crosslinked Polyethylene Liner, Neutral
		Longevity Liner, 7mm Offset
		Longevity Liner, Constrained
		Longevity Liner, Oblique
		Longevity Revision Shell Liner - 0° Neutral
		Maxera Cup
		Mayo 6 Deg Taper
		Mayo Stem
		Mayo Stem HA/TCP
		Metabloc Hip Stem cemented
		Metabloc Hip Stem uncemented
		Metasul CLS Spotorno Insert
		Metasul Femoral Head 12/14
		Metasul Femoral Head 8/10
		Metasul LDH Large Diameter Head
		Metasul LDH Large Diameter Head Adapter with 12/14 Cone
		Metasul LDH Large Diameter Head Adapter with 8/10 Cone
		Metasul Taper Liner
		Modular Bipolar Insert
		Modular Bipolar Shell
		Monoblock Acetabular Cup
		Monoblock Acetabular Cup - 10°
		Monoblock Acetabular Cup with Screw holes - 0°
		Monoblock Acetabular Cup with Screw holes - 10°
		Morscher Press-Fit cup, flat
		Morscher Press-Fit cup, Metasul
		Morscher Press-Fit cup, Sulene-PE
		MS-30 Distal Centralizer
		MS-30 Stem Lateral
		MS-30 Stem Standard
		Mueller Full Profile Cup
		Natural Hip (12/14 Taper) - Calcar Replacement
		Natural Hip (12/14 Taper) - CoCr NonPorous Offset
		Natural Hip (12/14 Taper) - Porous Collared Offset
		Natural Hip (12/14 Taper) - Porous Collared with Short neck
		Natural-Hip Stem, CoCr Slim
		Natural-Hip Stem, Non-Porous CoCr
		Natural-Hip Stem, Porous Collared

Manufacturer	Category	Implant Overlay Template Family
		Natural-Hip Stem, Porous Collarless
		Natural-Hip Stem, Porous Collarless Slim
		Natural-Hip stem, Porous with Collar Slim
		Nitrided CoCr Femoral Head
		Nitrided Femoral Head 6°
		Optan Stem Cemented
		Optan Stem Uncemented
		Original M.E. Mueller Acetabular Roof Reinforcement Ring
		Original M.E. Mueller CDH Stem D28, Protasul-S30
		Original M.E. Mueller Low Profile Cup, Durasul
		Original M.E. Mueller Low Profile Cup, Metasul
		Original M.E. Mueller Low Profile Cup, Sulene-PE
		Original M.E. Mueller Straight Stem lateral, Protasul-10
		Original M.E. Mueller Straight Stem lateral, Protasul-S30
		Original M.E. Mueller Straight Stem standard, Protasul-10
		Original M.E. Mueller Straight Stem standard, Protasul-S30
		PE CLS Spotorno Cup Insert
		PF Revision Stem Lateral
		PF Revision Stem Standard
		PF Stem Lateral
		PF Stem Lateral HA
		PF Stem Standard
		PF Stem Standard HA
		Poly Liners Constrained
		Porous with Collar Slim
		Pressfit Cup, Alpha
		Pressfit Cup, Kappa
		Pressfit Cup, Kappa HA
		Pressfit Cup, Kappa HA removable Plugs
		Protasul-S30 Femoral Head 12/14
		Protasul-S30 Femoral Head 8/10
		Reinforcement Ring with Hook
		Reinforcement Ring with Hook and Cranial Flange
		Revitan Distal Part Curved
		Revitan Distal Part, Straight
		Revitan Proximal Part, Conical
		Revitan Proximal Part, Cylindrical
		SL Titanium Shell
		St. Nabor Shell
		Staffit Cup
		Standard Polyethylene Liner, 10° Elevated
		Standard Polyethylene Liner, 20° Elevated
		Standard Polyethylene Liner, 7mm Offset
		Standard Polyethylene Liner, Eccentric
		Standard Polyethylene Liner, Neutral
		Standard Polyethylene Liner, Oblique
		Sulene PE Alpha Insert
		Sulene PE Alpha Insert hooded
		Sulox Ceramic Femoral Head 12/14
		Sulox Ceramic Femoral Head 8/10
		TPP Threaded Bolt

Manufacturer	Category	Implant Overlay Template Family
		TPP Thrust Plate
		Trabecular Metal Acetabular Augment
		Trabecular Metal Acetabular Shell, Cluster-Holed
		Trabecular Metal Acetabular Shell, Multi-Holed
		Trabecular Metal Acetabular Shell, Non-Holed
		Trabecular Metal Primary Hip
		Trabecular Metal Primary Hip Extended
		Trabecular Metal Revision Shell
		Tribosul Femoral head 12/14
		Trilogy AB Cluster-Hole Shell
		Trilogy AB Cluster-Hole Shell Cup
		Trilogy AB HA/TCP Multi-Hole Shell
		Trilogy AB Multi-Hole Shell
		Trilogy Acetabular Cup
		Trilogy F/M Acetabular Shell
		Trilogy F/M Acetabular Shell, Cluster-Holed
		Trilogy F/M Acetabular Shell, Multi-Holed
		Trilogy F/M Acetabular Shell, Non-Holed
		Trilogy F/M Acetabular Shell, Spiked
		Trilogy HA/TCP Acetabular Shell, Cluster-Holed
		Trilogy HA/TCP Acetabular Shell, Multi-Holed
		Trilogy HA/TCP Acetabular Shell, Non-Holed
		Trilogy HA/TCP Acetabular Shell, Spiked
		Trilogy HA/TCP Acetabular Shell, Uni-Holed
		Trilogy IT HA/TCP, Cluster-Hole
		Trilogy IT HA/TCP, Multi-Hole
		Trilogy IT, Cluster-Hole
		Trilogy IT, Multi-Hole
		Trilogy Longevity Constrained Liner
		Uni-Holed
		Unipolar Head 12/14
		Versys 10" Revision Stem Bowed +10, Beaded FullCoat Calcar
		Versys 10" Revision Stem Bowed +20, Beaded FullCoat Calcar
		Versys 10" Revision Stem Bowed, Beaded FullCoat
		VerSys 10" Revision Stem Straight, Beaded FullCoat
		VerSys 6" Stem Straight, Beaded FullCoat LM Extra Extended Offset
		VerSys 6" Stem Straight, Beaded FullCoat Low Head Center
		Versys 7.5" Revision Stem Straight +10, Beaded FullCoat Calcar
		Versys 7.5" Revision Stem Straight +20, Beaded FullCoat Calcar
		VerSys 7.5" Revision Stem Straight, Beaded FullCoat Plus
		Versys 8" Revision Stem Bowed, Beaded FullCoat Plus
		VerSys Advocate Hip Stem Non V-Lign, Extended Offset
		VerSys Advocate Hip Stem Non V-Lign, Standard Offset
		VerSys Advocate Hip Stem V-Lign, Extended Offset
		VerSys Advocate Hip Stem V-Lign, Standard Offset
		VerSys Cemented Revision Calcar
		VerSys Cemented Stem Short Neck, 6° Taper
		VerSys Collared 6" Beaded FullCoat Plus LM Body, Ext Offset
		VerSys Collared 6" Beaded FullCoat Plus LM Body, Std Offset
		VerSys Collared 6" Beaded FullCoat Plus Std Body, Ext Offset
		VerSys Collared 6" Beaded FullCoat Plus Std Body, Std Offset

Manufacturer	Category	Implant Overlay Template Family
		VerSys Collared 6° Beaded FullCoat Stem
		VerSys Collared Beaded MidCoat - Extra Ext Offset
		VerSys Collared Beaded MidCoat - Low Head Center Ext Offset
		VerSys Collared Cemented Body Stem
		VerSys Collared Cemented LD/Fx Stem
		VerSys Collared Fiber Metal MidCoat LM Body, Std Offset HA/TCP coated
		VerSys Collared Fiber Metal MidCoat Std Body, Std Offset HA/TCP coated
		VerSys Collared Fiber Metal MidCoat Stem
		VerSys Collared Press-Fit LD/Fx
		VerSys CRC Modular Build-Up Blocks
		VerSys Enhanced Taper Stem
		VerSys Epoch FullCoat Beaded Low Head Center (LHC) Stem
		VerSys Epoch FullCoat Body LM, Extended Offset Stem
		VerSys Epoch FullCoat Body LM, Extra-Extended Offset Stem
		VerSys Epoch FullCoat Body LM, Standard Offset Stem
		VerSys Epoch FullCoat Reduced Neck Length {RNL} 125 Stem
		VerSys Epoch FullCoat Reduced Neck Length {RNL} 135 Stem
		VerSys Femoral Head 12/14
		VerSys Fiber Metal Mid Coat Collared LM HA/TCP Coated, 6° Taper
		VerSys Fiber Metal Mid Coat Collared STD HA/TCP Coated, 6° Taper
		VerSys Fiber Metal Mid Coat Collarless LM HA/TCP Coated, 6° Taper
		VerSys Fiber Metal Mid Coat Collarless RM HA/TCP Coated, 6° Taper
		VerSys Fiber Metal Mid Coat Collarless STD HA/TCP Coated, 6° Taper
		VerSys Fiber Metal MidCoat Collarless LM Standard Offset HA/TCP coated
		VerSys Fiber Metal MidCoat Collarless Standard Offset HA/TCP coated
		VerSys Fiber Metal Taper LM Body, Std Offset HA/TCP coated
		VerSys Fiber Metal Taper LM Body, Std Offset HA/TCP Coated, 6° Taper
		VerSys Fiber Metal Taper LM Body, Std Offset, 6° Taper
		VerSys Fiber Metal Taper Std Body, Std Offset HA/TCP coated
		VerSys Fiber Metal Taper Std Body, Std Offset HA/TCP Coated, 6° Taper
		VerSys Fiber Metal Taper Std Body, Std Offset, 6° Taper
		VerSys Heritage CDH
		VerSys Heritage Primary Std Body, Ext Offset
		VerSys Heritage Primary Std Body, Std Offset
		VerSys Heritage Revision
		VerSys LD/FX Cemented, 6° Taper
		VerSys LD/FX Press Fit, 6° Taper
		Virtec Lateral Stem
		Virtec Long Stem
		Virtec Standard Stem
		Vitamin E Highly Crosslinked Liner, Elevated
		Vitamin E Highly Crosslinked Liner, Neutral
		Wagner Cone Prosthesis 125°
		Wagner Cone Prosthesis 135°
		Wagner SL Revision Stem
		Wagner SL Revision Stem 145°
		Wagner STANDARD Cup
		Weber Cup Metasul
		Weber Cup Sulene-PE
		Weber Stem, curved
		Weber Stem, straight

Manufacturer	Category	Implant Overlay Template Family
		ZCA All-Poly 10° Inclined Face Cup - 22mm
		ZCA All-Poly 10° Inclined Face Cup - 26mm
		ZCA All-Poly 10° Inclined Face Cup - 28mm
		ZCA All-Poly 10° Inclined Face Cup - 32mm
		ZCA All-Poly 10° Neutral Inclined Face Flanged Cup - 22mm
		ZCA All-Poly 10° Neutral Inclined Face Flanged Cup - 26mm
		ZCA All-Poly 10° Neutral Inclined Face Flanged Cup - 28mm
		ZCA All-Poly Neutral Cup – 22mm
		ZCA All-Poly Neutral Cup – 26mm
		ZCA All-Poly Neutral Cup – 28mm
		ZCA All-Poly Neutral Cup – 32mm
		ZCA All-Poly Neutral Flanged Cup - 22mm
		ZCA All-Poly Neutral Flanged Cup - 26mm
		ZCA All-Poly Neutral Flanged Cup - 28mm
		ZCA All-Poly Snap In Cup with Cement Spacers
		ZCA Longevity Crosslinked Poly 10° Inclined Face Cup
		ZCA Longevity Crosslinked Poly 10° Neutral Inclined Face Flanged Cup
		ZCA Longevity Crosslinked Poly Neutral Cup
		ZCA Longevity Crosslinked Poly Neutral Flanged cup
		ZCA Longevity Crosslinked Poly Snap In Cup with Cement Spacers
		Zimmer M/L Taper Extended Offset
		Zimmer M/L Taper Extended Offset HA/TCP
		Zimmer M/L Taper Extended Offset, Reduced Neck Length
		Zimmer M/L Taper Extended Offset, Reduced Neck Length HA/TCP
		Zimmer M/L Taper Standard Offset
		Zimmer M/L Taper Standard Offset HA/TCP
		Zimmer M/L Taper Standard Offset, Reduced Neck Length
		Zimmer M/L Taper Standard Offset, Reduced Neck Length HA/TCP
		Zimmer M/L Taper with Kinectiv Modular Neck Technology
		Zimmer M/L Taper with Kinectiv Modular Neck Technology HA/TCP
		Zimmer MMC Cup
		Zimmer MMC Cup HA
		Zimtron Femoral Head 12/14
		Zimtron Femoral Head 6°
		Zintra Threaded Cup
		ZMR Calcar A Body
		ZMR Calcar AA Body
		ZMR Calcar B Body
		ZMR Calcar C Body
		ZMR Calcar D Body
		ZMR Calcar E Body
		ZMR Calcar F Body
		ZMR Cone A Body
		ZMR Cone AA Body
		ZMR Cone B Body
		ZMR Cone Body
		ZMR Cone C Body
		ZMR Cone D Body
		ZMR Cone E Body
		ZMR Cone F Body
		ZMR Porous Stem Bowed

Manufacturer	Category	Implant Overlay Template Family
		ZMR Porous Stem Straight
		ZMR Slotted Porous Stem Bowed
		ZMR Splined Stem Bowed
		ZMR Splined Stem Straight
		ZMR Spout A Body
		ZMR Spout AA Body
		ZMR Spout B Body
		ZMR Spout C Body
		ZMR Spout D Body
		ZMR Spout E Body
		ZMR Spout F Body
		ZMR Taper Stem Straight
		ZMR XL Body Ext Offset
		ZMR XL Body Std Offset
		ZMR XL Porous Stem Bowed
		ZMR XL Porous Stem Straight
		ZMR XL Taper Stem
	Total Knees	Distal Femoral Component
		Fluted Stem Extension Bowed Precoat, cemented
		Fluted Stem Extension Straight Precoat
		Fluted Stem Mobile Tibial
		Gender Solutions Natural-Knee Flex Congruent Sulene Tibial Articular Surface
		Gender Solutions Natural-Knee Flex Female Non-Porous Femoral
		Gender Solutions Natural-Knee Flex Female Porous Femoral
		Gender Solutions Natural-Knee Flex Male Non-Porous Femoral
		Gender Solutions Natural-Knee Flex Male Porous Femoral
		Gender Solutions Natural-Knee Flex Prolong Congruent Tibial Articular Surface
		Gender Solutions Natural-Knee Flex Prolong Ultracongruent Tibial Articular Surface
		Gender Solutions Natural-Knee Flex Ultracongruent Sulene Tibial Articular Surface
		Innex Femoral Augment
		Innex Gender CR Femoral Implant, cemented
		Innex Gender CR Femoral Implant, cementless
		Innex Gender UCOR Mobile Articular Surface
		Innex Patella, cemented
		Intercalary Segment
		M/G Unicompartmental Articular Surface
		MBK CR/CS Articular Surfaces
		MBK CR/CS Porous Femorals
		MBK Option Fluted Stem Tibial Plate
		MIS CR-Flex Mobile Tibial Component
		MIS Drop Down Stem Extension
		Mobile Bearing Knee Femoral
		Natural-Knee Flex Ultracongruent Rotating Sulene Articular Surface
		Natural-Knee II Cobalt Chrome Modular Tibial Baseplate
		Natural-Knee II Congruent Tibial Articular Surface
		Natural-Knee II CR Non-Porous Femoral
		Natural-Knee II CR Porous Femoral
		Natural-Knee II CRK Non-Porous Femoral
		Natural-Knee II Modular Tibial Spacer
		Natural-Knee II Non-Porous Pegless Stemmed Tibial Baseplate
		Natural-Knee II Non-Porous Rotating Tibial Baseplate

Manufacturer	Category	Implant Overlay Template Family
		Natural-Knee II Porous Resurfacing Tibial Baseplate
		Natural-Knee II Porous Rotating Tibial Baseplate
		Natural-Knee II Porous Stemmed Tibial Baseplate
		Natural-Knee II Porous Tibial Spacer
		Natural-Knee II PS Non-Porous Femoral
		Natural-Knee II RC Femur Nonporous Revision Adapter Spacer
		Natural-Knee II Rev/Cons Femur Nonporous Distal Spacer
		Natural-Knee II Rev/Cons Femur Nonporous Posterior Spacer
		Natural-Knee II Revision Femoral Stem
		Natural-Knee II Revision Femoral Stem Bowed
		Natural-Knee II Revision Non-Porous Femoral
		Natural-Knee II Revision Non-Porous Femoral Augment
		Natural-Knee II Revision Non-Porous Tibial Baseplate
		Natural-Knee II Revision Porous Femoral
		Natural-Knee II Revision Porous Tibial Baseplate
		Natural-Knee II Revision Stem
		Natural-Knee II Revision Tibial Stem
		Natural-Knee II Revision/Constrained Tibial Insert
		Natural-Knee II Rotating Platform Congruent Insert
		Natural-Knee II Rotating Platform Ultra-Congruent Insert
		NexGen All-Poly Patella
		NexGen Cemented Stem Extension
		NexGen CR All-Poly Tibials
		NexGen CR Option Femoral
		NexGen CR Porous Femoral
		NexGen CR Porous Pegged Tibial Plate
		NexGen CR Precoat Femoral
		NexGen CR Precoat Pegged Tibial Plate
		NexGen CR/CR-Flex/CRA Prolong Articular Surface
		NexGen CRA Precoat Femoral
		NexGen CR-Flex Mobile
		NexGen CR-Flex Mobile Prolong Articular Surface
		NexGen CR-Flex Option Femoral
		NexGen CR-Flex Porous Femoral
		NexGen Extra Long Curved Stem Extension
		NexGen Extra Long Stem Extension
		NexGen Fluted 3° Option Stem Tibial Plate
		NexGen Fluted 7° Option Stem Tibial Plate
		NexGen Fluted Stem Mobile Precoat Tibial Plate
		NexGen Fluted Stem Tibial
		NexGen Full Block Tibial Augment
		NexGen Gender Solutions CR-Flex Porous Femoral
		NexGen Gender Solutions CR-Flex Precoat Femoral
		NexGen Gender Solutions LPS-Flex Option Femoral
		NexGen Gender Solutions LPS-Flex Plus Option Femoral
		NexGen Gender Solutions LPS-Flex Porous Femoral
		NexGen HA/TCP Tibial Plate, Pegged, Porous
		NexGen HA/TCP Tibial Plate, Stemmed, Porous
		NexGen LCCK Articular Surface
		NexGen LCCK Option Femoral
		NexGen Legacy 3° Fluted Precoat Tibial Augment

Manufacturer	Category	Implant Overlay Template Family
		NexGen LPS All Poly Tibial
		NexGen LPS Articular Surface
		NexGen LPS Femoral, Option
		NexGen LPS Option Femorals
		NexGen LPS Porous Femoral
		NexGen LPS Precoat Femoral
		NexGen LPS/LPS-Flex TM Monoblock Tibia
		NexGen LPS-Flex Fixed Articular Surface Highly Crosslinked
		NexGen LPS-Flex Fixed Articular Surface Machined
		NexGen LPS-Flex Micro (Net-Shape Molded) Articular Surface
		NexGen LPS-Flex Mobile Articular Surface
		NexGen LPS-Flex Option Femoral
		NexGen LPS-Flex Plus Femoral Component
		NexGen LPS-Flex Porous Femoral
		NexGen LPS-Flex Precoat Femoral
		NexGen LPS-Flex Titanium Femoral
		NexGen MIS LPS Mobile Tibial Plate
		NexGen MIS Modular Precoat Tibial Plate
		NexGen MIS Precoat Stemmed Tibial Plate
		NexGen Offset Long Stem Extension
		NexGen Offset Stem Extensions
		NexGen Option Stemmed Tibial Plate
		NexGen Porous Stemmed Tibial Plate
		NexGen Precoat 22deg Third Wedge Tibial Augment
		NexGen Precoat 7deg Full Wedge Tibial Augment
		NexGen Precoat A/P Wedge Stemmed Tibial Plate
		NexGen Precoat Anterior Femoral Augments
		NexGen Precoat CRA/LCCK/RHK Distal Femoral Augment
		NexGen Precoat CRA/LCCK/RHK Distal Only Femoral Augment
		NexGen Precoat CRA/LCCK/RHK Posterior Femoral Augment
		NexGen Precoat Posterior Only Femoral Augments
		NexGen Precoat Stemmed Tibial Plate
		NexGen Primary Porous Patella
		NexGen Prolong CR Flex Articular Surface
		NexGen RHK Hinge Post Extension
		NexGen Rotating Hinge Knee Articular Surface
		NexGen Rotating Hinge Knee Femoral
		NexGen Rotating Hinge Knee Full Block Tibial Augment Full Block
		NexGen Rotating Hinge Knee Precoat Tibial Plate
		NexGen Sharp Fluted Stem Extension
		NexGen Straight Stem Extensions
		NexGen Tibial Augment Block_15mm
		NexGen Tibial Augment Block_20mm
		NexGen Trabecular Metal Tibial Tray
		NexGen UC-Flex Mobile Articular Surface
		NexGen UC-Flex Mobile Prolong Articular Surface
		NextGen Precoat Tibial Augment
		Non-Stemmable Keel
		Persona CR Articular Surface
		Persona CR Narrow Femoral Component
		Persona CR Standard Femoral Component

Manufacturer	Category	Implant Overlay Template Family
		Persona PS Articular Surface
		Persona PS Narrow Femoral Component
		Persona PS Standard Femoral Component
		Persona Stemmed Cemented Tibia
		Persona Stemmed TM Tibia
		Persona Tapered Stem Extension
		Persona UC Articular Surface
		Segment with Male/Female Taper
		Segment with Male/Male Taper
		Segmental Articular Surface
		Segmental Basic Prox Femur 38mm Offset
		Segmental Smooth Collar
		Segmental TM Collar
		Segmental TM Femoral Diaphyseal Cone
		Segmental TM Femoral Metaphyseal Cone
		Segmental TM Proximal Femoral Curved Washer Kit
		Segmental TM Proximal Femoral Flat Washer Kit
		Segmental TM Proximal Femoral Offset
		Segmental TM Proximal Tibia
		Segmental TM Tibial Cone, Large
		Segmental TM Tibial Cone, Medium
		Segmental TM Tibial Tissue Washer Kit
		Stemmable Keel
		Tibial Augmentation Half Wedge
		Variable Stiffness Stem Extension Bowed Precoat
		Variable Stiffness Stem Extension Straight Precoat
		Zimmer Trabecular Metal Femoral Cone Augment
		Zimmer Trabecular Metal Tibial Cone
	Total Shoulders	Anatomical Shoulder Combined - Humeral Head Adaptor
		Anatomical Shoulder Fracture - Humeral Head with Base-Plate
		Anatomical Shoulder Fracture - Humeral Long Stem
		Anatomical Shoulder Fracture - Humeral Stem
		Anatomical Shoulder Inverse - Glenoid Fixation
		Anatomical Shoulder Inverse - Glenoid Head
		Anatomical Shoulder Inverse - Humeral cup
		Anatomical Shoulder Inverse - Inlay
		Anatomical Shoulder Inverse - Screw
		Anatomical Shoulder Removable Head - Glenoid with Keel
		Anatomical Shoulder Removable Head - Glenoid with Pegs
		Anatomical Shoulder Removable Head - Humeral Head
		Anatomical Shoulder Removable Head - Humeral Stem cemented
		Anatomical Shoulder Removable Head - Humeral Stem cementless
		Anatomical Shoulder Removable Head - Humeral Stem revision
		Bigliani/Flatow - Humeral Std. Head
		Bigliani/Flatow - Modular Humeral Stem
		Bigliani/Flatow - Modular Offset Head
		Bigliani/Flatow - Monoblock Humeral Stem
		Bigliani/Flatow Glenoid
		Bigliani/Flatow Humeral Offset Head
		Bigliani/Flatow Modular Stem Micro Implant
		Durom Shoulder - Cup

Manufacturer	Category	Implant Overlay Template Family
		Glenoid with Keel
		Glenoid with Keel Articular Surface
		Glenoid with Pegs
		Glenoid with Pegs Articular Surface
		Sidus Stem-Free Shoulder Humeral Anchor
		Sidus Stem-Free Shoulder Humeral Head
		TM Glenoid Articular Surface
		TM Reverse Base Plate
		Trabecular Metal Acetabular Restrictor
		Zimmer Trabecular Metal Reverse Shoulder - Dual Taper Insert
		Zimmer Trabecular Metal Reverse Shoulder - Glenosphere
		Zimmer Trabecular Metal Reverse Shoulder - Retentive Poly Liner 12° (65° inclination angle)
		Zimmer Trabecular Metal Reverse Shoulder - Spacer
		Zimmer Trabecular Metal Reverse Shoulder - Standard Poly Liner 7° (60° inclination angle)
		Zimmer Trabecular Metal Reverse Shoulder - TM Reverse Stem
		Zimmer Trabecular Metal Shoulder - Humeral Stem, 48°
		Zimmer Trabecular Metal Shoulder - TM Stem 42°
		Trauma
	4.0mm NCB Cannulated Screw self-drilling self-tapping, sterile	
	4.0mm NCB Screw self-tapping, non-sterile	
	4.0mm NCB Screw self-tapping, sterile	
	4.5mm NCB Cancellous Screw self-tapping, non-sterile	
	4.5mm NCB Cancellous Screw self-tapping, sterile	
	5.0mm NCB Cancellous Screw self-tapping thread length 32mm, non-sterile	
	5.0mm NCB Cancellous Screw self-tapping thread length 32mm, sterile	
	5.0mm NCB Screw self-tapping, non-sterile	
	5.0mm NCB Screw self-tapping, sterile	
	6.5mm Cervical Screw cannulated non-sterile	
	6.5mm Cervical Screw cannulated sterile	
	Antegrade GT Femoral Nail Ti, sterile	
	Antegrade PF Femoral Nail Universal Ti, sterile	
	Cable-Ready Bone Plates	
	Cable-Ready Integral Crimps	
	Cancellous Screw Ø 6.0, Titanium, Fixed Angle 3.5mm Hex Head, non-sterile	
	Cancellous Screw Ø 6.0, Titanium, Fixed Angle 3.5mm Hex Head, sterile	
	Cancellous Screw Ø 6.0, Titanium, Partially Threaded 3.5mm Hex Head, non-sterile	
Cancellous Screw Ø 6.0, Titanium, Partially Threaded 3.5mm Hex Head, sterile		
Cephalomedullary Femoral Nail - Asia Ti-6AL-4V, sterile		
Cephalomedullary Femoral Nail - Long Stainless Steel, sterile		
Cephalomedullary Femoral Nail - Long Ti-6AL-4V, sterile		
Cephalomedullary Femoral Nail - Short Stainless Steel, sterile		
Cephalomedullary Femoral Nail - Short Ti-6AL-4V, sterile		
Cephalomedullary Nail Cap SST, non-sterile		
Cephalomedullary Nail Cap SST, sterile		
Cephalomedullary Nail Cap Ti-6AL-4V, non-sterile		
Cephalomedullary Nail Cap Ti-6AL-4V, sterile		
Cortical Screw Ø 4.0, Titanium, Fixed Angle 3.5mm Hex Head, non-sterile		
Cortical Screw Ø 4.0, Titanium, Fixed Angle 3.5mm Hex Head, sterile		
Cortical Screw Ø 4.0, Titanium, Partially Threaded 3.5mm Hex Head, non-sterile		
Cortical Screw Ø 4.0, Titanium, Partially Threaded 3.5mm Hex Head, sterile		
Cortical Screw Ø 5.0, Titanium, Fixed Angle 3.5mm Hex Head, sterile		

Manufacturer	Category	Implant Overlay Template Family
		Femoral Nail Cap Titanium, non-sterile
		Femoral Nail Cap Titanium, sterile
		ITST Intramedullary System Femoral Nail
		ITST Intramedullary System Femoral Nail Short
		ITST Intramedullary System Lag Screws
		ITST Nail Caps
		Lag Screw Stainless Steel, non-sterile
		Lag Screw Stainless Steel, sterile
		Lag Screw Ti-6AL-4V, non-sterile
		Lag Screw Ti-6AL-4V, sterile
		Lag Screw, non-sterile
		M/DN Intramedullary System Femoral Nail
		M/DN Intramedullary System Humeral Nail
		M/DN Intramedullary System Retrograde Femoral Nail
		M/DN Intramedullary System Tibial Nail
		M/DN Screws
		NCB Curved Femur Shaft Plate, non-sterile
		NCB Curved Femur Shaft Plate, sterile
		NCB Distal Femoral Plate, non-sterile
		NCB Distal Femoral Plate, sterile
		NCB Periprosthetic Distal Femur plate, non-sterile
		NCB Periprosthetic Distal Femur plate, sterile
		NCB Periprosthetic Proximal Femur plate, non-sterile
		NCB Periprosthetic Proximal Femur plate, sterile
		NCB Periprosthetic Trochanter Plate, sterile
		NCB Proximal Humeral Plate with Polyaxial Locking Screws, sterile
		NCB Proximal Humeral Plate, sterile
		NCB Proximal Tibial Plate 2L, non-sterile
		NCB Proximal Tibial Plate 2L, sterile
		NCB Proximal Tibial Plate 3L, non-sterile
		NCB Proximal Tibial Plate 3L, sterile
		NCB Screw Ø 4.0, deep thread, non-sterile
		NCB Screw Ø 4.0, deep thread, sterile
		NCB Tuberculum Minus Plate, non-sterile
		NCB Tuberculum Minus Plate, sterile
		NCB Unicortica Ø 5.0, non-sterile
		NCB Unicortica Ø 5.0, sterile
		Periarticular Plating System – 4.5mm Distal Medial Femoral plate
		Periarticular Plating System – Distal Dorsal Delta Radial plate
		Periarticular Plating System – Distal Dorsal T Radial Plate
		Periarticular Plating System – Distal Lateral Femoral Plate
		Periarticular Plating System – Distal Lateral Fibular Plate
		Periarticular Plating System – Distal Lateral Tibial Plate
		Periarticular Plating System – Distal Medial Humeral plate
		Periarticular Plating System – Distal Medial Tibial Plate
		Periarticular Plating System – Distal Posterior Lateral Fibular Plate
		Periarticular Plating System – Distal Volar Radial Radius Plate
		Periarticular Plating System – Distal Volar Ulnar Radius Plate
		Periarticular Plating System – Lateral Calcaneal Plate
		Periarticular Plating System – Posterior Lateral Humeral Plate
		Periarticular Plating System – Proximal Dorsal Ulnar Plate

Manufacturer	Category	Implant Overlay Template Family
		Periarticular Plating System – Proximal Lateral Humeral Plate
		Periarticular Plating System – Proximal Lateral Tibial Plate
		Periarticular Plating System – Proximal Medial Tibial Plate
		Retrograde Femoral Nail Cap Titanium, sterile
		Retrograde Femoral Nail Universal Titanium, sterile
		Retrograde Short Femoral Universal Titanium, sterile
		Sirus Intramedullary Nail for Femur D10.3 left, sterile
		Sirus Intramedullary Nail for Femur D10.3 right, sterile
		Sirus Intramedullary Nail for Femur D11 left, sterile
		Sirus Intramedullary Nail for Femur D11 right, sterile
		Sirus Intramedullary Nail for Femur D12 left, sterile
		Sirus Intramedullary Nail for Femur D12 right, sterile
		Sirus Intramedullary Nail for Femur D13 left, sterile
		Sirus Intramedullary Nail for Femur D13 right, sterile
		Sirus Intramedullary Nail for Femur D9.3 left, sterile
		Sirus Intramedullary Nail for Femur D9.3 right, sterile
		Sirus Intramedullary Nail for Proximal Humerus, sterile
		Sirus Intramedullary Nail for Tibia D10 cannulated, sterile
		Sirus Intramedullary Nail for Tibia D11 cannulated, sterile
		Sirus Intramedullary Nail for Tibia D12 cannulated, sterile
		Sirus Intramedullary Nail for Tibia D8.3 solid, sterile
		Sirus Intramedullary Nail for Tibia D9.3 cannulated, sterile
		Ti Compression Screw, non-sterile
		Tibia Nail Universal Ti, sterile
		Tibial Nail Cap Titanium, non-sterile
		Tibial Nail Cap Titanium, sterile
		Ti-Versa FxII Compression Tube/Plate with collar, Long Type, non-sterile
		Ti-Versa FxII Compression Tube/Plate with collar, Short Type, non-sterile
		Ti-Versa FxII Supracondylar Compression Tube/Plate, non-sterile
	Unicondylar Knees	Allegretto Femoral Implant Distally Reinforced
		Allegretto Femoral Implant Standard
		Allegretto Tibial Plate Metalback
		Allegretto Tibial Plate All-Polyethylene
		Gender Solutions PFJ Precoat
		Innex CR Fixed Articular Surface
		Innex CR Mobile Articular Surface
		Innex Femoral Spacer
		Innex Fixed Tibial Plate, cemented
		Innex Mobile Tibial Plate, cementless
		Innex Mobile Tibial Plate, cemented
		Innex SC Femoral Implant, cemented
		Innex SC Fixed Articular Surface
		Innex Standard CR Femoral Implant, cemented
		Innex Standard CR Femoral Implant, cementless
		Innex Stem Extension, Angle 2°
		Innex Stem Extension, Conical
		Innex Stem Extension, Offset
		Innex Stem Extension, Straight
		Innex Tibial Augment
		Innex Tibial Spacer
		Innex UC Fixed Articular Surface

Manufacturer	Category	Implant Overlay Template Family
		Innex UCOR Mobile Articular Surface
		M/G Unicompartmental All-Poly Tibial
		M/G Unicompartmental Precoat Femoral
		M/G Unicompartmental Tibial Plate
		Natural Knee II Unicompartmental Nonporous Tibial Baseplate
		Natural Knee II Unicompartmental Porous Tibial Baseplate
		Natural-Knee II Porous Unicompartmental SinterLock Femoral
		Natural-Knee II Unicompartmental Congruent Tibial Insert
		Natural-Knee II Unicompartmental Nonporous Femoral
		NexGen CR/CR-Flex TM Monoblock Tibia
		Persona Partial Knee - Medial Articular Surface (new)
		Persona Partial Knee - Medial Cemented Femoral Component (new)
		Persona Partial Knee - Medial Cemented Tibia Component (new)
		Zimmer High Flex Unicompartmental Articular Surface
		Zimmer High Flex Unicompartmental Tibial Plate
		Zimmer UniCompartmental High Flex All Poly Tibials - Lateral
		Zimmer UniCompartmental High Flex All Poly Tibials - Medial
		Zimmer UniCompartmental Precoat Femoral